

Grant

Special Offer

PCV Series Centrifuge/Vortexer Combination

Thomas No. 1228J30



List ~~\$700~~

Promo Price **\$500**

You Save **\$200!**

- Centrifugation speed: fixed at 2800 rpm/50Hz and 3500 rpm/60Hz
- Continuous operation or short spin
- Choice of interchangeable rotors for different microtube sizes/combinations and for 0.2 ml strips
- 2-Year manufacturer warranty

Cost-effective, fixed-speed, combined mini centrifuge/vortex mixer for combined or independent centrifuge and mixing applications of microtubes and 0.2 ml microtube strips in low volume applications.

Tubes are loaded into the rotor for simultaneous spinning, then removed for individual mixing in the vortex cup located at the top of the central shaft. In spin-mix-spin applications, the Combi-spin™ can be used with very low reagent volumes representing an overall saving in time, labor and material.

While Supplies Last!

Thomas
Scientific

We Believe You Are Important
— How Can We Help?

800.345.2100 • 800.345.5232 fax
PO Box 99 • Swedesboro, NJ 08085
www.thomassci.com • value@thomassci.com

