

FORENSICS

CONSTELLATIONS TAKES
CF RECONSTRUCTION
INTERDISCIPLINARY
EVIDENCE SCENE
EARTHQUAKE
FLUIDS
INVESTIGATIVE
DACTYLOSCOPY
DISPENSING
ELECTRANTS
ENTRANCE
S
CIRCUMSTANCES
COMPUTATIONAL
GAIN
FIELD
ARCHAEOLOGICAL
FINGERPRINTING
VEIN
EFFECTIVENESS
PAINT DICTIONARIES
OBJECTS QUESTIONED
PROFESSIONALLY GLASSES
CRIMINOLOGISTS POSSIBLE
PRODUCTS
COUNTERPARTS
CE
ISTS
UCT

Thomas Scientific Forensics 2015 Catalog Supplement



ARROWHEAD FORENSICS
DISCOVER WHAT YOU'VE BEEN MISSING

Universal Swab Boxes

These universal swab boxes are made of 18pt, virgin fiber, solid bleached sulfate. This stock helps prevent cross contamination and is the sturdiest swab box available. These 6.5" boxes are pre-cut with 2 breathing holes and designed to fit all swabs. Contains no biodegradable components to jeopardize the evidence. The breathable boxes are ideal for cotton tipped swabs, Cap-shure™ swabs and our Bio-Swab that have been used to collect body fluids. Ideal for swabs in sexual assault cases. Simply squeeze the box and the box pops up to prevent contamination.



Thomas No.	Mfr. No.	Description	Pack
1182X33	A-7034	Plain Universal Swab Box, 6.5"	100

Adhesive Evidence Labels

These labels include all information needed to clearly mark the evidence contained in the packaging. Labels are available with or without the chain of custody.



Thomas No.	Mfr. No.	Description	Pack
1182X24	A-3706	Adhesive Evidence Labels - Evidence, 4 x 6"	100
1182X25	A-3700	Adhesive Evidence Labels - Evidence, Chain of Custody, 4 x 6"	100

Zenith Magnetic Applicator

This highly magnetic wand is capable of carrying more than three times the powder of standard wands and features an innovative burnishing head to improve the quality of crime scene latent lifts. Features a spring-loaded ejection for easy release of powder. The pentagon shaped head ensures that the tip never touches the ground and prevents the applicator from rolling. The head is also removable so you can decontaminate between scenes.



Thomas No.	Mfr. No.	Description	Each
1182X37	A-2917	Zenith Magnetic Applicator	1

Sharps Evidence Collection Tubes

Exposure to blood-borne diseases and protection of personnel is a serious problem. These clear plastic tubes will preserve the evidence and protect evidence collection personnel. All three sizes include a custom foam insert for immobilization, evidence strips for sealing, and a biohazard label to clearly mark the evidence as a possible hazard.



Thomas No.	Mfr. No.	Description	Pack
1182X03	A-1907	Syringe Tube, 7 x 1"	12
1182X04	A-1905	Syringe Tube, 9 x 1"	12
1182X05	A-1918	Knife Tube, 18 x 3"	6

SealGuard™ Split Back Evidence Tape

Arrowhead SealGuard™ Split Back Evidence Tape features an acetate base to prevent shrinking and swelling in extreme temperatures. 3M adhesive produces a superior, tamper-proof seal, making it impossible to remove the tape without destroying the tape. Dual serrated edges make it easy to tear to desired length. SealGuard™ is backed with a 70 lb. split back liner to prevent curling, allowing quick, easy, controlled application.



Thomas No.	Mfr. No.	Description	Each
1182X38	A-3101	Red Solid "Evidence" Imprint, 1-3/8" x 108'	1
1182X39	A-3001	Red/White Stripe "Evidence" Imprint, 1-3/8" x 108'	1

DNA Caution Labels

With the increase of DNA as an element in crime-solving, preservation of DNA on collected evidence is becoming more critical. Make sure collected evidence is protected with our exclusive DNA Caution Labels. Available in two sizes to fit any item collected, the labels are a great addition to any crime scene evidence collection kit. Our DNA Caution Labels are printed in fluorescent green and black on self-adhesive stock with peel-away backing, making them easily visible.



Thomas No.	Mfr. No.	Description	Roll
1182X26	A-3731	DNA Caution Label, 1 x 1.5"	100
1182X27	A-3732	DNA Caution Label, 2.5 x 2"	100

Arrowhead Forensics (Continued)

"Patrol" Standard Fingerprint Kit

Many departments are now providing basic fingerprint kits for patrol cars and requiring those officers to lift fingerprints. Arrowhead's Patrol Kits contain everything needed for these officer's to do this.

Each Kit Contains:

- 1 - 2 oz. Standard Black Powder
- 1 - Fiber Brush
- 1 - Roll 1-1/2" Clear Lifting Tape
- 50 - White Back Cards



Thomas No.	Mfr. No.	Description	Each
1182X36	A-2219W	"Patrol" Standard Fingerprint Kit	1

Wide Diameter Evidence Collection Tubes

These clear plastic tubes will preserve evidence and protect collection personnel. The 2" diameter is ideal for glass pipes. Each container includes a foam insert, security seal, and biohazard label.

Thomas No.	Mfr. No.	Description	Pack
1182X06	A-1917	Evidence Collection Tube, 7 x 2"	12



Latent Fingerprint Powders

You can be confident that the powder you get will be consistent in quality every time. All of our powders are precision milled and sifted to assure quality.



Thomas No.	Mfr. No.	Description	Each
1182X11	A-2310W	Black Powder, 2 oz.	1
1182X12	A-2312W	Black Powder, 16 oz.	1
1182X13	A-2320W	Dual Use Powder, 2 oz.	1
1182X14	A-2322W	Dual Use Powder, 16 oz.	1
1182X15	A-2324W	Dual Use Powder, 128 oz.	1
1182X16	A-2340W	White Powder, 2 oz.	1
1182X17	A-2344W	White Powder, 128 oz.	1

Magnetic Latent Fingerprint Powders

All of our magnetic powders are made with the same care and attention to detail as our standard powders. This powder must be applied with a magnetic applicator that acts as a "brush" by creating a ball of powder. This ball is the only thing that touches the active surface.



Thomas No.	Mfr. No.	Description	Each
1182X18	A-2412W	Black Powder, 16 oz.	1
1182X19	A-2422W	Dual Use Powder, 16 oz.	1
1182X20	A-2442W	White Powder, 16 oz.	1

Master Fingerprint Kit

Each Kit Contains:

- 2 - 2 oz. Black Powder
- 1 - 4 oz. Black Magnetic Powder
- 1 - 4 oz. Dual Use Magnetic Powder
- 1 - 4 oz. Fluorescent Blue Magnetic Powder
- 1 - 4 oz. Fluorescent Orange Magnetic Powder
- 1 - 4 oz. Fluorescent Yellow Magnetic Powder
- 1 - 4 oz. Fluorescent Pink Magnetic Powder
- 1 - roll 2" Frosted Lifting Tape Powder
- 1 - roll 2" Clear Lifting Tape
- 1 - Mini Scalpel Cutter
- 1 - Bottle Rough Lift
- 2 - Fiber Brushes
- 1 - Feather Brush
- 1 - Magnetic Applicator
- 50 - White Lifting Cards
- 50 - Black Lifting Cards
- 1 - UV Lamp



Thomas No.	Mfr. No.	Description	Each
1182X23	A-2640B	Master Fingerprint Kit	1

Master Trace Evidence Collection Kit

Each Kit Contains:

- 10 - Straight Edge Razor Blades
- 15 - Antimicrobial Hand Wipes
- 40 - Swab Boxes
- 10 - Tongue Depressors
- 10 - Cuticle Sticks
- 5 - Combs
- 5 - 3.25 x 2 x .875" Evidence Slide Box
- 12 - Plastic Disposable Tweezers
- 5 - Round Blade Disposable Scalpels
- 5 - Pointed Blade Disposable Scalpels
- 50 - 6 x 9" Kraft Envelopes
- 50 - 2.75 x 4.25" Glassine Envelopes
- 10 - Transfer Pipettes
- 3 - 50ml Polypropylene Collection Tubes
- 40 - 6" Sterile Cotton Tip Swabs 2/Pack
- 2 - 5ml Sterile Water Ampules
- 2 - 2 Safety Goggles
- 2 - Black Sharpies
- 2 - Black Ink Pens
- 24 - 4 x 4" Clear J Lifts
- 10 - 4 x 6" 2ml Reclosable Evidence Bags
- 10 - 6 x 8" 4ml Reclosable Evidence Bags
- 50 - 2.875" x 5.25" Kraft Envelopes



Thomas No.	Mfr. No.	Description	Each
1182X34	A-8864	Master Trace Evidence Collection Kit	1

Backing Cards

After you have taken the time and effort to locate and lift a print at a crime scene, you must have the proper backing card to both preserve and document the print. Our backing cards are available in snow white or jet black on front. The backs of the cards have an area for all the information needed: case number, date, time, and diagram for the location of the print.



Thomas No.	Mfr. No.	Description	Pack
1182X45	A-2731	Backing Cards, White, 3 x 5"	100
1182X46	A-2732	Backing Cards, Black, 3 x 5"	100

Arrowhead Forensics (Continued)

Master Evidence Collection Kit

This Master Evidence Kit can be used in large agencies for multiple CSIs or one single kit for a small department. Packaged in a sturdy Pelican carry case with custom dividers and a pocketed pouch system. Built strong to withstand the demands from scene to scene.



To view entire kit contents, please visit: thomassci.com Keyword: **1182X34**

Thomas No.	Mfr. No.	Description	Each
1182X34	A-8868	Master Evidence Collection Kit	1

Master Crime Scene Tac Pak Kit

The all new Master Crime Scene Tac-Pak Kit features the products needed to accomplish basic and some advanced levels of crime scene processing. Secure the premises with barrier tape, mark evidence and from there move into the collection process the scene requires. Protect valuable GSR and trace evidence with our signature Hand Preservation Bags, seal evidence on the spot with our Seal-Guard Tamper Resistant Evidence Tape and collect DNA samples with the Bio-Swab. The kit will give you the ability to collect latent fingerprints, trace evidence, impression evidence, digital evidence and biologicals. Proper evidence packaging for each of these collected items is provided. The black tactical bag comes with a shoulder strap for easy transport.



To view entire kit contents, please visit: thomassci.com Keyword: **1182X44**

Thomas No.	Mfr. No.	Description	Each
1182X44	A-8866	Master Crime Scene Tac Pak Kit	1

Bio-Swab™

Bio-Swab™ is a pre-filled, individually wrapped, gamma radiated, sterile applicator that will deliver the exact amount of water to collect a DNA sample. This packaging system eliminates the need for carrying sterile water and will further protect from cross contamination. The Bio-Swab™ can be transported to the lab in air dry boxes or coin envelopes. Simply hold the applicator with the color ring end pointing upwards, gently snap the tip with the color coded end, and allow the sterile water to fill the applicator. Validation material available upon request.



Thomas No.	Mfr. No.	Description	Pack
1182X31	A-7013	Bio-Swab Sterile Applicator	100
1182X32	A-7014BP	Bio-Swab Sterile Applicator - Bulk Pack	500

Single Use Latent Print Fiber Brush

Our fiberglass single use latent print brush is designed for crime scene investigators concerned with cross contamination from scene to scene. Hand made with the same fibers as our standard brush with a sturdy plastic handle. Each brush is individually packaged in a wrapper.



Thomas No.	Mfr. No.	Description	Each
1182X22	A-2498	Single Use Latent Print Fiber Brush	1

Hand Preservation Bags

This breathable bag can be used to cover the hands of a suspected shooter, resulting in the preservation of valuable GSR evidence. This bag can also be used to prevent contamination on the hands of the deceased. Made of Tyvek® material with a drawstring and new string-lock attachment for securing the bag in place without tying.



Thomas No.	Mfr. No.	Description	Pair
1182X40	A-4830	7 x 12-1/2"	1
1182X41	A-4830P	7 x 12-1/2"	50

Tie-Down Weapons Storage Boxes

These boxes are made of 200 lb. corrugated cardboard and are designed to keep evidence from shifting in their boxes. They will help you protect prints on the evidence by not allowing the prints to be smudged in transit. They are also excellent for storing large numbers of items and keeping your property room organized. Features tie down holes for a secure fastening of weapons.



Thomas No.	Mfr. No.	Description	Pack
1182W75	A-1060	Knife Boxes & Ties, 16 x 3 x 2" (25 Boxes/50 Ties)	25
1182W76	A-1050	Long Gun Boxes & Ties, 50-1/2 x 6-3/4 x 3" (25 Boxes/50 Ties)	25
1182W77	A-1070	Pistol Boxes & Ties, 12 x 7-1/2 x 2" (25 Boxes/50 Ties)	25

Cyclone Fiber Brush

Our best Fiber brush available, the Cyclone has been precision milled to spin, swirl and evenly deliver powder over the desired surface without damaging the print. 6.5" long brush consists of thousands of specially formulated fibers to ensure softness and maximum coverage without damaging prints. The brush will hold the correct amount of powder and distribute it evenly, resulting in fewer trips to the powder jar. Our Fiber Brushes are made to last, and will grow better with age. Packaged in a reusable plastic tube which opens on both ends for easy storage and removal. Use one brush for black powder, one for dual use powder and one for white powder.



Thomas No.	Mfr. No.	Description	Each
1182X21	A-2488	Cyclone Fiber Brush, 6.5" long	1

3M Fingerprint Lifting Tape

Unquestionably the finest quality tape available for lifting prints. 3M tape is pressure wound, heat treated, and fish-eye free.



Thomas No.	Mfr. No.	Description	Each
1182X07	A-2101	Fingerprint Lifting Tape, Clear, 1 x 360"	1
1182X08	A-2102	Fingerprint Lifting Tape, Clear, 1-1/2 x 360"	1
1182X09	A-2103	Fingerprint Lifting Tape, Clear, 2 x 360"	1
1182X10	A-2104	Fingerprint Lifting Tape, Clear, 4 x 360"	1

Arrowhead Forensics (Continued)

The Monster™ Evidence Bags

These bags are ideal for collecting and storing oversized evidence such as bed linens, pillows, rugs, clothing, shoes, and other saturated or moist evidence. These Monster™ Bags may easily be sealed with evidence tape. The Monster™ Bag is made with 2-ply, 50 lb., gusseted Kraft construction paper that is extremely strong and durable. It is freestanding and resistant to tearing.



Thomas No.	Mfr. No.	Description	Pack
1182X42	A-1034	16-1/2 x 34 x 11-1/2" (50 Bags & 50 Labels)	50
1182X43	A-1034 NFR	16-1/2 x 34 x 11-1/2" (50 Bags)	50

Paper Evidence Bags – Printed

Simply place your evidence inside the bag and seal with evidence seal or tape. You can log your information along with a chain of possession on the outside of the bag. Each printed bag is clearly marked and includes areas for required evidence submission information band chain of custody. Paper bags are easily sealed with our Seal-Guard™ Evidence Tape or Evidence Box Sealing Tape. The convenient sizes allow for the packaging of a variety of evidence items including clothing, shoes, and drugs. Because the paper allows for the evidence to dry, they are the best choice for packaging DNA sensitive evidence. Available in three sizes.



Thomas No.	Mfr. No.	Description	Pack
1182W66	A-1012	Paper Evidence Bags, Printed, 7 x 4 x 13-1/2" (12 lb.)	100
1182W67	A-1012C	Paper Evidence Bags, Printed, 7 x 4 x 13-1/2" (12 lb.)	500
1182W68	A-1025	Paper Evidence Bags, Printed, 8 x 5-1/4 x 15-1/2" (25 lb.)	100
1182W69	A-1025C	Paper Evidence Bags, Printed, 8 x 5-1/4 x 15-1/2" (25 lb.)	500
1182W70	A-1016	Paper Evidence Bags, Printed, 12 x 7 x 17" (57 lb.)	100
1182W71	A-1016C	Paper Evidence Bags, Printed, 12 x 7 x 17" (57 lb.)	500

Paper Evidence Bags – Unprinted

Sturdy Kraft paper bags in various sizes. These unprinted bags have space for the information you need to properly mark your evidence. Unprinted bags are available for use with separate labels. We recommend sealing these bags with SealGuard™ evidence tapes.



Thomas No.	Mfr. No.	Description	Pack
1182W72	A-1012P	Paper Evidence Bags, Unprinted, 7-1/16 x 4-1/2 x 13-3/4" (12 lb.)	100
1182W73	A-1025P	Paper Evidence Bags, Unprinted, 8-1/4 x 5-15/16 x 15-7/8" (25 lb.)	100
1182W74	A-1016P	Paper Evidence Bags, Unprinted, 12 x 7 x 17" (57 lb.)	100

SureSeal™ Tamper Indicating Currency Bags

These 3 mil thick currency/coin bags maintain the same features and integrity as the personal property bags but are made exclusively for holding seized currency and coins. Arrowhead's currency bag has a preprinted label to clearly identify the money it contains, as well as why the money was collected.



Thomas No.	Mfr. No.	Description	Pack
1182W86	A-1163-M	Tamper Indicating Currency Bag, 6 x 9"	100
1182W87	A-1163-MC	Tamper Indicating Currency Bag, 6 x 9"	500
1182W88	A-1193-MC	Tamper Indicating Currency Bag, 9 x 12"	500

SureSeal™ Tamper Indicating Evidence Bags w/Control Numbers

NEW!

These 3 mil Sure-Seal Tamper Indicating bags require no prepackaging of evidence prior to use. Once the adhesive closure strip is exposed and sealed it is permanent and impossible to open without destroying the seal and showing obvious signs of tampering. Includes white write-on area, tear and puncture resistance durability, removable submission receipt, and the convenience of clear, see-through material. Available in four sizes.



Thomas No.	Mfr. No.	Description	Pack
1182W78	A-1143	Tamper Indicating Evidence Bag, w/Control Number, 4 x 9"	100
1182W79	A-1143BP	Tamper Indicating Evidence Bag, w/Control Number, 4 x 9"	500
1182W80	A-1163	Tamper Indicating Evidence Bag, w/Control Number, 6 x 9"	100
1182W81	A-1163BP	Tamper Indicating Evidence Bag, w/Control Number, 6 x 9"	500
1182W82	A-1193	Tamper Indicating Evidence Bag, w/Control Number, 9 x 12"	100
1182W83	A-1193BP	Tamper Indicating Evidence Bag, w/Control Number, 9 x 12"	500
1182W84	A-1123	Tamper Indicating Evidence Bag, w/Control Number, 12 x 16"	100
1182W85	A-1123BP	Tamper Indicating Evidence Bag, w/Control Number, 12 x 16"	500

VaporLok Heavy Duty 5 mil Nylon Barrier Heat Seal Bags

These Pouches are ideal for packaging and storing chemicals, controlled substances, paraphernalia, and volatile materials. Barrier Pouches are useful in a wide variety of applications for evidence collection and feature a 3/8" 3 sided seal. These 5 mil heat sealable nylon pouches are constructed of a combination of films with ideal barrier properties. This barrier eliminates cross-contamination while protecting individuals from exposure to dangerous gasses "leaking" through. These pouches are approved and used by many federal, state, and local narcotics agencies. Pouches may be sealed with any quality Heat Sealer.



Thomas No.	Mfr. No.	Description	Pack
1182W89	A-1174	Heavy Duty Nylon Barrier Heat Seal Bags, 5 mil, 4 x 6"	100
1182W90	A-1174C	Heavy Duty Nylon Barrier Heat Seal Bags, 5 mil, 4 x 6"	1000
1182W91	A-1176	Heavy Duty Nylon Barrier Heat Seal Bags, 5 mil, 6-1/2 x 8"	100
1182W92	A-1176C	Heavy Duty Nylon Barrier Heat Seal Bags, 5 mil, 6-1/2 x 8"	1000
1182W93	A-1178	Heavy Duty Nylon Barrier Heat Seal Bags, 5 mil, 8 x 10"	100
1182W94	A-1178C	Heavy Duty Nylon Barrier Heat Seal Bags, 5 mil, 8 x 10"	1000

For additional information and package sizes, please visit: thomasci.com Keyword: 1182W89

SealGuard™ Biohazard Square Labels

The SealGuard™ Biohazard Square labels contain the international Biohazard Symbol and are printed on fluorescent orange stock to insure visibility and are provided on a roll for easy dispensing.



Thomas No.	Mfr. No.	Description	Roll
1182X28	A-3738	Biohazard Square Label, 1 x 1"	250
1182X29	A-3735	Biohazard Square Label, 1-1/2" x 1-1/2"	250

SealGuard™ Large Biohazard Labels

The following labels contain the international Biohazard Symbol and are printed on fluorescent orange stock to insure visibility and are provided on a roll for easy dispensing.

Thomas No.	Mfr. No.	Description	Pack
1182X30	A-3720	Large Biohazard Label, 4 x 5"	100





Connecticut Analytical Corporation

High Speed Latent Fingerprint Kit

FastPrint is a simple, rapid fingerprinting kit that can develop latent fingerprints from any porous substance or surface in only 5 minutes.



Thomas No.	Mfr. No.	Description	Each
1180R02	1005-HHDNFK1	High Speed Latent Fingerprint Kit	1
1180R03	1005-HHDNFK1R	High Speed Latent Fingerprint Material Refill Kit	1

Dried Blood Spot or Specimen (DBS) Spotting Cards

Chemical free, general purpose high speed analytical-fast metabolite DNA and pharmacokinetic blood spotting cards.

Blue cards offer 4 positions to spot blood.

Red cards offer 12 positions to spot blood with a convenient hydrophobic protective flap.



Thomas No.	Mfr. No.	Description	Pack
1180R02	P1375-CR01250	Dried Blood Spot or Specimen (DBS) Blue Spotting Cards	100
1180R03	P1375-CB04100	Dried Blood Spot or Specimen (DBS) Red Spotting Cards	50

Graphene Oxide Dried Blood Spot Extraction Buffer Solution

Unleash "The Power of Graphene". Graflex DBS Extraction Buffer can triple the extraction concentration of multiple analytes from whole blood samples within minutes. Specifically designed for dried blood spot analysis for clinical and R&D use only. For best results, we recommend using our HSA-FMDP spotting cards.



Thomas No.	Mfr. No.	Description	Each
C855E32	1005-81001	Graphene Oxide Dried Blood Spot Extraction Buffer Solution	1

Graphene Oxide DNA Extraction Buffer Solution

Unleash "The Power of Graphene". Graflex DNA Extraction Buffer is designed to detect DNA and nucleic acids for sample analysis in clinical and R&D laboratories. This patent-pending solution has been optimized to specifically target and bind to nucleic acids with high affinity and without causing interference, nor includes detergents of any kind that can degrade DNA. Scientifically proven to obtain a 3 fold increase in DNA extraction versus other elution and extraction buffers.



Thomas No.	Mfr. No.	Description	Each
C855E33	1005-81002	Graphene Oxide DNA Extraction Buffer Solution	1

Graphene Oxide Proteomics Extraction Buffer Solution

Unleash "The Power of Graphene". Graflex Proteomics Extraction Buffer can significantly increase the extraction concentration of peptides, proteins and protein targets within minutes. Specifically designed for pharmacokinetic analysis in clinical and R&D laboratories. Our patent-pending Graphene Oxide solution has scientifically shown more target proteins and peptides were extracted versus other extraction methods. Graflex Proteomics Extraction Buffer is perfect for those doing drug development research and clinical dosaging studies.



Thomas No.	Mfr. No.	Description	Each
C855E34	1005-81003	Graphene Oxide Proteomics Extraction Buffer Solution	1

Graphene Oxide Chemistry Extraction Buffer Solution

Unleash "The Power of Graphene". Graflex Chemistry Extraction Buffer can significantly increase the extraction concentration of peptides, proteins and protein targets within minutes. Specifically designed for pharmacokinetic analysis in clinical and R&D laboratories. Our patent-pending Graphene Oxide solution has scientifically shown more target proteins and peptides were extracted versus other extraction methods.



Thomas No.	Mfr. No.	Description	Each
C855E35	1005-81004	Graphene Oxide Chemistry Extraction Buffer Solution	1

Graphene Oxide Forensic Extraction Buffer Solution

Unleash "The Power of Graphene". Like Graflex DNA solution, Graflex Forensics Extraction Buffer triples the extraction concentration of analytes from whole blood and other biological samples in just minutes. This patent-pending Graphene Oxide solution is ideal for recovery of multiple biomarkers and biological analytes from specimen samples such as blood, urine, tears, saliva and stool.



Thomas No.	Mfr. No.	Description	Each
C855E36	1005-81005	Graphene Oxide Forensic Extraction Buffer Solution	1



Micro-Chem Plus™

- Hospital grade / Insures complete broad spectrum disinfecting
- Dual-quat / Residual bactericidal action
- Kills AIDS virus / Reduces danger of contamination from hard surfaces
- Biodegradable, phosphate free / Environmentally responsible
- Multiple applications / Variety of uses for both institutional and industrial areas




Thomas No.	Mfr. No.	Description	Case
C849T34	0255-29	Micro-Chem Plus, 4 x 1 gal	4



Integra® Miltex® Finger Ring Cutter

For safe and efficient ring removal.



Thomas No.	Mfr. No.	Description	Each
1177K08	33-140	Finger Ring Cutter, with circular saw blade, 6-1/4" (160mm)	1



Dissecting Set – Pro

Stainless steel, consisting of 20 instruments for accurate dissection. Supplied in zipper leather case.

Each Kit Contains:

- Scissors - 2 pcs.
- Forceps - 2 pcs.
- Needles - 4 pcs.
- Scalpel - 1 pc.
- B.P. Handle w Blade - 1 pc.
- Razor - 1 pc.
- Ruler 6" - 1 pc.
- Magnifying Glass - 1 pc.
- Glass Slide - 1 pc.
- Chain and Hook - 1 pc.
- Cover Slips - 1 pc.
- Blow Pipe - 1 pc.
- Needle Wooden Handle - 1 pc.
- Dropper Glass - 1 pc.
- Brush - 1 pc.



Thomas No.	Mfr. No.	Description	Each
1177L67	BI0148	Dissecting Set – 20 instruments	1

Dissecting Set – Popular

Stainless steel, consisting of 14 instruments for accurate dissection. Complete in zipper leather case.

Each Kit Contains:

- Scissors - 2 pcs.
- Forceps - 2 pcs.
- Needles - 4 pcs.
- Scalpel - 1 pc.
- B.P. Handle w Blade - 1 pc.
- Razor - 1 pc.
- Magnifying Glass - 1 pc.
- Chain and Hook - 1 pc.
- Blow Pipe - 1 pc.



Thomas No.	Mfr. No.	Description	Each
1177L68	BI0150	Dissecting Set – 14 instruments	1

Dissecting Set – Beginner

Stainless steel, consisting of 7 instruments for accurate dissection. Complete in zipper leather case.

Each Kit Contains:

- Scissors - 2 pcs.
- Forceps - 2 pcs.
- Seeker - 1 pc.
- B.P. Handle w/Blade - 1 pc.
- Needle Wooden Handle - 1 pc.



Thomas No.	Mfr. No.	Description	Each
1177L69	BI0152	Dissecting Set – 7 instruments	1

Dissecting Set – Economy

Stainless steel, consisting of 7 instruments for accurate dissection. Complete in zipper leather case.

Each Kit Contains:

- Scissors - 1 pcs.
- Forceps - 1 pcs.
- Needles - 2 pcs.
- B.P. Handle w/Blade - 1 pc.
- Dropper Glass - 1 pc.
- Ruler 6" - 1 pc.



Thomas No.	Mfr. No.	Description	Each
1177L70	BI0152ECO	Dissecting Set – Economy, 7 instruments	1

Dissecting Set – Student

Stainless steel, consisting of 9 instruments for accurate dissection. Complete in zipper leather case.

Each Kit Contains:

- Scissors - 1 pcs.
- Forceps - 1 pcs.
- Needle Wooden Handle - 2 pcs.
- Scalpel - 1 pc.
- Scalpel Handle w/Blade - 1 pc.
- Dropper Glass - 1 pc.



Thomas No.	Mfr. No.	Description	Each
1177L71	BI0149	Dissecting Set – Student, 9 instruments	1

Dissecting Set – Senior

Stainless steel, consisting of 9 instruments for accurate dissection. Complete in zipper leather case.

Each Kit Contains:

- Scissors - 2 pcs.
- Forceps - 2 pcs.
- Needle - 2 pcs.
- Scalpel - 1 pc.
- Scalpel Handle w/Blade - 1 pc.
- Ruler 6" - 1 pc.



Thomas No.	Mfr. No.	Description	Each
1177L72	BI0151	Dissecting Set – Senior, 9 instruments	1

Dissecting Set – Advanced

Stainless steel, consisting of 12 instruments for accurate dissection. Complete in zipper leather case.

Each Kit Contains:

- Scissors - 2 pcs.
- Forceps - 2 pcs.
- Needle - 2 pcs.
- Artery Forceps - 1 pc.
- Tissue Forceps - 1 pc.
- Seeker - 1 pc.
- Scalpel Handle w/ Blade - 1 pc.
- Bone Cutter - 1 pc.
- Dental Probe & Seeker - 1 pc.



Thomas No.	Mfr. No.	Description	Each
1177L73	BI0151ADV	Dissecting Set – Advanced, 12 instruments	1

UNITED SCIENTIFIC SUPPLIES, INC.

Ink Chromatography and Forensics STEM Kit

Correlated to the Next Generation Science Standards

Use this kit to introduce your students to ink chromatography and its application in solving forensics cases. Students gain experience in creating ink databases, identifying 'ink tags' to 'date' various inks and then use them (as professional document examiners do) in solving two mysteries involving document inks – an apparent suicide note and a forgery case. Includes materials for three detailed activities (with 'model' and 'inquiry' investigational paths – a total of four investigational activities), which can accommodate up to 40 students. Also contains a DVD with PDF Teacher and Student Guides and other digital content. Covered concepts include: data analysis, chromatographs, inks and their composition-mixtures, compounds, mixtures, molecules, solubility, affinity, polar substances, Rf values, capillary action, light energy, fluorescence, and hypothesis testing.



Thomas No.	Mfr. No.	Description	Each
1201U30	AISCRKIT	Ink Chromatography and Forensics STEM Kit	1

Ultraviolet Light

This high-quality 6", 4W handheld blacklight is ideal for forensics or illuminating fluorescent and black-light reactive items. Also includes a built-in LED flashlight. Requires 4 AA batteries (not included).



Thomas No.	Mfr. No.	Description	Each
1177J56	AISUULT1	Ultraviolet Light	1



Sterile Forensic Applicators

Overall Length: 5.125" (130.18 mm)

Handle: Polypropylene

Tip: Large Round

Tip Material: Foam

Sterile

Applications/Characteristics:

- Crime scene investigation
- Buccal cell collection
- Specimen collection
- Features large swab surface for optimum use
- Individually wrapped Puritan Pur-Wraps brand



Thomas No.	Mfr. No.	Description	Case
1215C53	25-1805 1PF RND FDNA	Large Foam Tip, Polypropylene Shaft, Sterile	500
1207Y91	25-3306 U TT FDNA	DNA Free, Ultra Flocked Swab	500
1213Q28	25-806 2PD TT	6" Poly, Sterile, Transport Tube, 24/30/2	1440
1215C52	25-806 2WC TT FDNA	Cotton Tip, Wood Shaft, Transport Tube	500
1213Q32	25-807 1PC EC	7" Cotton Swab, Cap-shure, 10/50/1	500

Kimberly-Clark

KLEENGUARD® A60 Bloodborne Pathogen and Chemical Splash Protection Sleeve Protectors

Protection against bloodborne pathogens in forensic or other environments where exposure to blood is possible.

Features:

- Elastic top
- Strong for extended wear
- Three-layer fabric construction features a middle layer of microporous film that allows heat and sweat vapor to escape while protecting against dry particulates and liquid
- Passes ASTM F1670/1671 testing for penetration of blood, body fluids and bloodborne pathogens
- 21 inches in length



Thomas No.	Mfr. No.	Description	Case
1184A28	36871	KLEENGUARD A60 Bloodborne Pathogen and Chemical Splash Protection Sleeve Protector, 21", Blue	200

KLEENGUARD® A60 Bloodborne Pathogen and Chemical Splash Protection Coveralls

Protection against bloodborne pathogens in forensic or other environments where exposure to blood is possible.

Features:

- REFLEX® design
- Three-layer fabric construction features a middle layer of microporous film that allows heat and sweat vapor to escape while protecting against dry particulates and liquid
- Passes ASTM F1670/1671 testing for penetration of blood, body fluids and bloodborne pathogens
- Seamless front provides more protection in primary exposure areas
- Taped storm flap
- Outer layer provides protection of film barrier against tearing and abrasion



Thomas No.	Mfr. No.	Description	Case
1232G24	45002	KLEENGUARD A60 Bloodborne Pathogen and Chemical Splash Protection Coveralls, Denim Blue, M	24
1232G25	45003	KLEENGUARD A60 Bloodborne Pathogen and Chemical Splash Protection Coveralls, Denim Blue, L	24
1232G26	45004	KLEENGUARD A60 Bloodborne Pathogen and Chemical Splash Protection Coveralls, Denim Blue, XL	24
1184A29	45005	KLEENGUARD A60 Bloodborne Pathogen and Chemical Splash Protection Coveralls, Denim Blue, 2XL	24
1184A30	45006	KLEENGUARD A60 Bloodborne Pathogen and Chemical Splash Protection Coveralls, Denim Blue, 3XL	20
1232G27	45007	KLEENGUARD A60 Bloodborne Pathogen and Chemical Splash Protection Coveralls, Denim Blue, 4XL	20

SIGMA-ALDRICH

AUTHORIZED DISTRIBUTOR

Titanium Tweezers

- Lightweight (high strength to weight ratio)
- Zero magnetic retentively
- Temperature resistant to 800°F
- Chemical and corrosion resistant
- Non-solder adhesive
- Bio-compatible (maintain cell integrity, no inflammatory response)



Thomas No.	Mfr. No.	Description	Each
1179V61	Z680192-1EA	Titanium tweezers, Style 00	1
1179V62	Z680257-1EA	Titanium tweezers, Style 6	1
1179V63	Z680265-1EA	Titanium tweezers, Style 7	1

Thomas Scientific Pipet Tips

Tips for truly accurate sample delivery

Generate reproducible results in minimal time, with minimal fuss. Our convenient, quality-tested pipette tips consistently deliver precise sample measurements in single and multi-channel pipettors:

- Virtually eliminate sample retention with our low retention polymer pipette tips
- Get reliable results with pipette tips that are certified free of RNase, DNase and endotoxins
- Confirm your draw visually through our graduated sides



Pre-Sterile Racked Pipet Tips

Thomas No.	Description	Pack	Case
1182T07	10µL, Pre-Sterile, Racked (10 Racks of 96 per Pack)	960	4800
1182T08	10XL, Pre-Sterile, Racked (10 Racks of 96 per Pack)	960	4800
1182T09	200µL, Pre-Sterile, Racked (10 Racks of 96 per Pack)	960	4800
1182T10	1000µL, Pre-Sterile, Racked (8 Racks of 96 per Pack)	768	3840
1182T11	1250µL, Pre-Sterile, Racked (8 Racks of 96 per Pack)	768	3840

Non-Sterile Racked Pipet Tips

Thomas No.	Description	Pack	Case
1182T12	10µL, Non-Sterile, Racked (10 Racks of 96 per Pack)	960	4800
1182T13	10XL, Non-Sterile, Racked (10 Racks of 96 per Pack)	960	4800
1182T14	200µL, Non-Sterile, Racked (10 Racks of 96 per Pack)	960	4800
1182T15	1000µL, Non-Sterile, Racked (8 Racks of 96 per Pack)	768	3840
1182T16	1250µL, Non-Sterile, Racked (8 Racks of 96 per Pack)	768	3840

Pre-Sterile Racked Filter Pipet Tips

Thomas No.	Description	Pack	Case
1182T02	10µL, Pre-Sterile, Racked (10 Racks of 96 per Pack)	960	4800
1182T03	10XL, Pre-Sterile, Racked (10 Racks of 96 per Pack)	960	4800
1182T04	200µL, Pre-Sterile, Racked (10 Racks of 96 per Pack)	960	4800
1182T05	1000µL, Pre-Sterile, Racked (8 Racks of 96 per Pack)	768	3840
1182T06	1250µL, Pre-Sterile, Racked (8 Racks of 96 per Pack)	768	3840

Bulk Pipet Tips

Thomas No.	Description	Pack	Case
1181F63	200µL, Graduated, Low Retention, Bulk (1000 tips/pack)	1000	10000
1182T28	10µL, Bulk (1000 tips/pack)	1000	20000
1182T29	200µL, Bulk (1000 tips/pack)	1000	10000
1182T30	200µL, Yellow, Bulk (1000 tips/pack)	1000	10000
1182T31	1000µL, Bulk (1000 tips/pack)	1000	4000
1182T32	1000µL, Blue, Bulk (1000 tips/pack)	1000	4000

Non-Sterile Low Retention Reload Pipet Tips

Thomas No.	Description	Pack	Case
1182T22	10µL, Non-Sterile, Low Retention, Reloads	960	4800
1182T23	10XL, Non-Sterile, Low Retention, Reloads	960	4800
1182T24	200µL, Non-Sterile, Low Retention, Reloads	960	4800
1182T25	1000µL, Non-Sterile, Low Retention, Reloads	768	3840
1182T26	1250µL, Non-Sterile, Low Retention, Reloads	768	3840

Pre-Sterile Low Retention Racked Pipet Tips

Thomas No.	Description	Pack	Case
1182T17	10µL, Pre-Sterile, Low Retention, Racked (10 Racks of 96 per Pack)	960	4800
1182T18	10XL, Pre-Sterile, Low Retention, Racked (10 Racks of 96 per Pack)	960	4800
1182T19	200µL, Pre-Sterile, Low Retention, Racked (10 Racks of 96 per Pack)	960	4800
1182T20	1000µL, Pre-Sterile, Low Retention, Racked (8 Racks of 96 per Pack)	768	3840
1182T21	1250µL, Pre-Sterile, Low Retention, Racked (8 Racks of 96 per Pack)	768	3840

Pre-Sterile Low Retention Racked Filter Pipet Tips

Thomas No.	Description	Pack	Case
1182R95	10µL, Pre-Sterile, Low Retention, Racked (10 Racks of 96 per Pack)	960	4800
1182R96	10XL, Pre-Sterile, Low Retention, Racked (10 Racks of 96 per Pack)	960	4800
1182R97	20µL, Pre-Sterile, Low Retention, Racked (10 Racks of 96 per Pack)	960	4800
1182R98	100µL, Pre-Sterile, Low Retention, Racked (10 Racks of 96 per Pack)	960	4800
1182R99	200µL, Pre-Sterile, Low Retention, Racked (10 Racks of 96 per Pack)	960	4800
1182T00	1000µL, Pre-Sterile, Low Retention, Racked (8 Racks of 96 per Pack)	768	3840
1182T01	1250µL, Pre-Sterile, Low Retention, Racked (8 Racks of 96 per Pack)	768	3840

Sterile Petri Dishes

Generic Sterile Petri Dishes are available in stackable and slippable forms

- Optically clear
- Vented lids for air circulation
- Radiation sterilized

These Stackable Petri Dishes are ideal for general laboratory use while Slippable Petri Dishes are designed for automated filling systems. Vented lids permit air circulation during incubation to reduce condensation build-up. Petri dishes are optically clear for easy viewing.

Petri Dishes are available in the following sizes: 60 x 15 mm, 100 x 15 mm (Stackable and Slippable), and 150 x 15 mm



Thomas No.	Description	Case
1188N80	Petri Dish, 100 x 15mm, Slippable, Sterile, 500	500
1188N81	Petri Dish, 100 x 15mm, Stackable, Sterile, Bulk, 25/500	500
1188N82	Petri Dish, 60 x 15mm, Stackable, Sterile, 25/500	500
1188N83	Petri Dish, 150 x 15mm, Slippable, Sterile, Bulk, 10/100	100