

MATERIAL SAFETY DATA SHEET



IBI SCIENTIFIC
9861 KAPP COURT
PEOSTA, IA 52068

EMERGENCY TELEPHONE NUMBER:
3E COMPANY (800) 451-8346
(001) 760-602-8703

TELEPHONE NUMBER FOR INFORMATION:
(800) 253-4942 OR (563) 690-0484
E-MAIL: INFO@IBISCI.COM

DATE PREPARED: 9-4-98

DATE REVISED: 1-13-10

SECTION 1. ----- CHEMICAL IDENTIFICATION -----

CATALOG#: IB01110
NAME: 10X TRIS-GLYCINE BUFFER CONCENTRATE
CAS#: N/A BLEND

SECTION 2. ----- HAZARDS IDENTIFICATION -----

IRRITANT

IRRITATING TO EYES, RESPIRATORY SYSTEM AND SKIN.
IN CASE OF CONTACT WITH EYES, RINSE IMMEDIATELY WITH PLENTY OF
WATER AND SEEK MEDICAL ADVICE.
WEAR SUITABLE PROTECTIVE CLOTHING.

HMIS: NFPA:
H: 1 H: N/A
F: 0 F: N/A
R: 0 R: N/A
PE: E SH: N/A

SECTION 3. ----- COMPOSITION/INFORMATION ON INGREDIENTS -----

CHEMICAL NAME:	CAS#:	RTECS #:
TRIS	77-86-1	TY2900000

FOR MORE INFORMATION SEE COMPLETE RTECS ENTRY.

SECTION 4. ----- FIRST-AID MEASURES -----

IN CASE OF CONTACT, IMMEDIATELY FLUSH EYES WITH COPIOUS AMOUNTS OF
WATER FOR AT LEAST 15 MINUTES.
IN CASE OF CONTACT, IMMEDIATELY WASH SKIN WITH SOAP AND COPIOUS
AMOUNTS OF WATER.
IF INHALED, REMOVE TO FRESH AIR. IF NOT BREATHING GIVE ARTIFICIAL
RESPIRATION. IF BREATHING IS DIFFICULT, GIVE OXYGEN.
IF SWALLOWED, WASH OUT MOUTH WITH WATER PROVIDED PERSON IS CONSCIOUS.
CALL A PHYSICIAN.
WASH CONTAMINATED CLOTHING BEFORE REUSE.

SECTION 5. ----- FIRE FIGHTING MEASURES -----

EXTINGUISHING MEDIA
WATER SPRAY.
CARBON DIOXIDE, DRY CHEMICAL POWDER OR APPROPRIATE FOAM.
SPECIAL FIREFIGHTING PROCEDURES
WEAR SELF-CONTAINED BREATHING APPARATUS AND PROTECTIVE CLOTHING TO
PREVENT CONTACT WITH SKIN AND EYES.

UNUSUAL FIRE AND EXPLOSIONS HAZARDS
EMITS TOXIC FUMES UNDER FIRE CONDITIONS.

SECTION 6. ----- ACCIDENTAL RELEASE MEASURES-----

WEAR RESPIRATOR, CHEMICAL SAFETY GOGGLES, RUBBER BOOTS AND HEAVY RUBBER GLOVES.
COVER WITH DRY LIME OR SODA ASH, PICK UP, KEEP IN A CLOSED CONTAINER AND HOLD FOR WASTE DISPOSAL.
VENTILATE AREA AND WASH SPILL SITE AFTER MATERIAL PICKUP IS COMPLETE.

SECTION 7. ----- HANDLING AND STORAGE-----

AVOID INHALATION.
AVOID CONTACT WITH EYES, SKIN AND CLOTHING.
AVOID PROLONGED OR REPEATED EXPOSURE.
WASH THOROUGHLY AFTER HANDLING.
KEEP TIGHTLY CLOSED.
STORE IN A COOL DRY PLACE.
STORE AT ROOM TEMPERATURE.

SECTION 8. ----- EXPOSURE CONTROLS/PERSONAL PROTECTION-----

CHEMICAL SAFETY GOGGLES.
COMPATIBLE CHEMICAL-RESISTANT GLOVES.
NIOSH/MSHA-APPROVED RESPIRATOR.
SAFETY SHOWER AND EYE BATH.
MECHANICAL EXHAUST RECOMMENDED.

SECTION 9. ----- PHYSICAL AND CHEMICAL PROPERTIES -----

APPEARANCE
LIQUID

SECTION 10. ----- STABILITY AND REACTIVITY -----

STABILITY
STABLE.
INCOMPATIBILITIES
STRONG OXIDIZING AGENTS
HAZARDOUS COMBUSTION OR DECOMPOSITION PRODUCTS
TOXIC FUMES OF:
CARBON MONOXIDE, CARBON DIOXIDE
NITROGEN OXIDES
AMMONIA
HAZARDOUS POLYMERIZATION
WILL NOT OCCUR.

SECTION 11. ----- TOXICOLOGICAL INFORMATION -----

ACUTE EFFECTS
MAY BE HARMFUL BY INHALATION, INGESTION, OR SKIN ABSORPTION.
CAUSES EYE AND SKIN IRRITATION.
MATERIAL IS IRRITATING TO MUCOUS MEMBRANES AND UPPER RESPIRATORY TRACT.
TO THE BEST OF OUR KNOWLEDGE, THE CHEMICAL, PHYSICAL, AND TOXICOLOGICAL PROPERTIES HAVE NOT BEEN THOROUGHLY INVESTIGATED.

SECTION 12. ----- ECOLOGICAL INFORMATION -----

DATA NOT YET AVAILABLE.

SECTION 13. ----- DISPOSAL CONSIDERATIONS -----

OBSERVE ALL FEDERAL, STATE AND LOCAL ENVIRONMENTAL REGULATIONS.

SECTION 14. ----- TRANSPORT INFORMATION -----

NOT REGULATED

SECTION 15. ----- REGULATORY INFORMATION -----

EUROPEAN INFORMATION

R: 36/37/38

IRRITATING TO EYES, RESPIRATORY SYSTEM AND SKIN.



Irritant (Xi)

S: 36/37/39

WEAR SUITABLE PROTECTIVE CLOTHING, GLOVES AND EYE/FACE PROTECTION.

SECTION 16. ----- OTHER INFORMATION -----

THE ABOVE INFORMATION IS BELIEVED TO BE CORRECT BUT DOES NOT PURPORT TO BE ALL INCLUSIVE AND SHALL BE USED ONLY AS A GUIDE. IBI SCIENTIFIC SHALL NOT BE HELD LIABLE FOR ANY DAMAGE RESULTING FROM HANDLING OR FROM CONTACT WITH THE ABOVE PRODUCT. SEE REVERSE SIDE OF INVOICE OR PACKING SLIP FOR ADDITIONAL TERMS AND CONDITIONS OF SALE.