

**CALL-IN
SPECIAL ONLY!
800.345.2100**

2013 EXCESS INVENTORY SALES PROGRAM

PFL® Headband Facemask



- Positive Facial Lock® Technology
- Magic Arch® Supports
- Breathe Easier & Safer
- 3-Ply Mask
- Delta P: 2.2 mm H₂O/cm²
- Latex Particle Challenge: 97.8%

Part No.
1204T96
(300/case)

\$70.00

PFL® Ear Loop Facemask



- Positive Facial Lock® Technology
- Magic Arch® Supports
- Breathe Easier & Safer
- 3-Ply Mask
- Delta P: 2.2 mm H₂O/cm²
- Latex Particle Challenge: 97.8%

Part No.
1204T97
(300/case)

\$98.00

PFL® Headband Facemask



- Positive Facial Lock® Technology
- No Magic Arch® Support
- Hypoallergenic
- 3-Ply Mask
- Delta P: 2.2 mm H₂O/cm²
- Latex Particle Challenge: 97.8%

Part No.
1204T98
(300/case)

\$68.00

PFL® Ear Loop Facemask



- Positive Facial Lock® Technology
- No Magic Arch® Support
- Hypoallergenic
- 3-Ply Mask
- Delta P: 2.2 mm H₂O/cm²
- Latex Particle Challenge: 97.8%

Part No.
1204T99
(300/case)

\$95.00

AlphaAir® Ear Loop Facemask



- Outstanding Particle Filtration
- Plastic Coated Nosepiece
- Soft Knit Ear Loops
- 3-Ply Mask
- Delta P: 2.3 mm H₂O/cm²
- Latex Particle Challenge: 97.5%

Part No.
1204U03
(500/case)

\$19.00

GenPro® Isolation Gown



- Moderate Barrier Protection
- Strong, Lightweight
- For Non-Hazardous Environments
- Full Coverage
- Ties at Neck & Waist

Part No.
1204U04
(50/case)

\$28.00

AlphaAir® Ear Loop Facemask



- Outstanding Particle Filtration
- Aluminum Nosepiece
- 3-Ply Mask
- Delta P: 4.4 mm H₂O/cm²
- Latex Particle Challenge: 93.5%
- Available in Blue, Pink and Yellow

Part No.
1204U02
(500/case)

\$19.00

Part No.
1204U00
1204U01
(500/case)

\$27.00

• **PFL® (Positive Facial Lock®)** - this feature utilizes the exclusive **Twist Seal®** technology that allows the user to create a custom fit each and every time

• **Magic Arch®** - this feature holds the mask away from the nose and mouth to create a larger breathing chamber to enhance user comfort

• **Non-Fluid Resistant Facemasks**

* While Supplies Last. All Sales Final, Non-Returnable. To receive this special pricing, please call 800.345.2100 to place your order.