


Notebooks

Laboratory

- Tamper-evident, section-sewn binding
- Acid-free, archival-safe paper
- Imprinted with 1/4" scientific grid or scientific ruling
- Available in 120 or 216 page formats

Thomas Laboratory and Engineering Notebooks feature the highest quality construction and are designed to protect your intellectual property. Notebooks feature a durable, hardbound, black imitation leather cover with special coating to protect lab work and prevent tampering. Notebooks include Guidelines for proper Laboratory Notebook documentation of your research, Table of Contents pages for recording page contents and dates and signature blocks for recording names, company, address, phone and email of owner. Special fields include: Title, Project, Book Number, "Continued from Book No." and "Continued to Book No." Laboratory Notebooks have a highly flexible round binding that allows the book to lay completely flat when open. Laboratory notebooks come in either grid or ruled format and pages measure 8-7/8" x 11-1/4".



Thomas No.	Description	Qty	Case
1236E02	Laboratory Notebook, Ruled, 120 Pg.	1	10
1236E04	Laboratory Notebook, Ruled, 216 Pg.	1	10
1236E06	Laboratory Notebook, Grid, 120 Pg.	1	10
1236E08	Laboratory Notebook, Grid, 216 Pg.	1	10

Continued on reverse side

Engineering

Engineering Notebooks are similar to Laboratory versions but possess square binding with section-sewn pages, ensuring that pages will not detach in rugged lab and field environments. Each notebook is imprinted with 1/4" scientific grids and include guidelines for proper documentation of your notebook. Notebook pages measure 8" x 10".



Thomas No.	Description	Qty	Case
1236E10	Engineering Notebook, Grid, 120 Pg.	1	10
1236E12	Engineering Notebook, Grid, 216 Pg.	1	10

Laboratory and Engineering Notebook Holders

- Side mount, under shelf and free standing styles
- Mount instantly with 3M® double stick tape

Side mount book holders can be mounted on the side of instruments, hoods or file cabinets. Under shelf model can be installed in a vertical or horizontal position under a bench or shelves. The free-standing holders offer sturdy placement for books on bench, desk or instrument.



Thomas No.	Description	Each
7084M88	Book Holder, Side Mount, Large, 6" x 2" x 6"	1
7084M94	Book Holder, Free Standing, 6" x 3" x 6"	1
7084M96	Book Holder, Free Standing, 6" x 4" x 6"	1
7084M98	Book Holder, Free Standing, 6" x 5" x 6"	1