

Safety SPOTLIGHT

Thomas
Scientific

We Believe You Are Important
— How Can We Help?

www.thomassci.com

Think **Thomas** for All
of Your **Safety** Needs!



Phone **800.345.2100**
Fax **800.345.5232**

Summer Edition 2017

Use **Promo Code 0217** to receive the featured pricing. Valid until 12/31/17.

Thomas
Scientific

Workrite
FR

Workrite
FR CP

THE RIGHT FIT FOR EVERY LAB

Innovations in Laboratory Safety



Workrite FR offers different styles of lab coats to meet various needs. All are flame-resistant (FR), but there are other key differences. Newer styles are made with a lighter-weight fabric, making them more breathable and comfortable. They also have knit cuffs, which were added based on feedback from the University of California. These cuffs keep foreign material out and prevent accidental spills.

The latest innovation is the Workrite FR/CP Lab Coat, which combines the FR properties of DuPont Nomex® IIIA with the chemical-splash protection (CP) of Westex ShieldCXP™. The first of its kind, this lab coat simultaneously protects against short-term thermal hazards and inadvertent liquid chemical splashes while also providing superior comfort and flexibility.

Chemicals that ShieldCXP is generally effective against include aqueous chemicals and mixtures, and polar organic solvents. For more information, refer to workritefrcp.com

Featuring fabrics from:

**DuPont™
Nomex®**  **Westex
ShieldCXP™**

© Workrite Uniform Company, Inc.
Workrite Uniform is an affiliate of Williamson-Dickie Manufacturing Company
Nomex® is a registered trademark of DuPont. | ShieldCXP™ is a trademark of Westex, Inc.

LAB COAT
Thomas Nos.
1169L19 – 1169L28



**LAB COAT
W/KNIT CUFFS**
Thomas Nos.
1169L39 – 1169L48



**FR/CP
LAB COAT**
Thomas Nos.
1169L29 – 1169L38



Flame-resistant (FR) protection	✓	✓	✓
Styles for men and women	✓	✓	✓
Concealed snap front closure	✓	✓	✓
Pass-through access to inner clothing		✓	✓
Knit cuffs for improved comfort and protection		✓	✓
Lightweight 4.5-ounce DuPont Nomex® IIIA fabric		✓	✓
Chemical-splash protection (CP) with Westex ShieldCXP™			✓

1.800.345.2100 | www.thomassci.com | value@thomassci.com



Think Thomas for All of Your Safety Needs!

Thomas
Scientific



New Pig

PIG® Spill Kit in Bucket

- Polyethylene bucket with easy-open, screw-on lid can be used to ship or temporarily store used absorbents
- PIG BLUE Socks stop spreading spills; PIG Mat absorbs quickly
- Prepacked container provides easy access to contents
- Large, color-coded label makes identification fast and accurate
- Tamperproof seals help prevent pilfering of spill response supplies
- UN Rating 1H2/Y30/S for Packing Groups II and III
- Use as refill for wall-mount holder (sold separately)
- Only PIG Spill Kits feature PIG Absorbents proven for rapid response.
- For information on custom spill kits, just call 1-800-HOT-HOGS (468-4647)



Thomas No.	Description	Each	List
1205U97	PIG® Spill Kit in Bucket	1	73.33
1205U98	PIG® Oil-Only Spill Kit in Bucket	1	85.99

Thomas
Scientific

Universal Spill Kit

- Neutralizes acids, bases organics, hydrocarbons, water reactive liquids, poisons, air reactive liquids, oxidizers, reducing agents and a host of other types of chemicals
- Color-coded and designed to change colors
- Non-toxic
- Safe to handle
- Easy-to-store
- Easy-to-use
- Rust inhibitor prevents degradation of hazmat shipping and storage containers



Ampho-Mag™ is a rapid response powder that safely adsorbs and treats any spill and produces less heat and spattering than any other neutralizer available. A gas neutralizer and adsorbent, it prevents pressure build-up in containers.

Thomas No.	Description	Each	List	Sale
8238N05	Universal Spill Kit	1	193.67	141.92



Looking for signage?

Thomas Scientific offers a wide variety of signs from National Marker (NMC) to keep your workplace safe.

Custom signage also available!

Ask your Thomas Scientific Sales Representative!

Visit thomassci.com
Search National Marker



Sure-Grip® EX Under-counter Corrosives/Acid Steel Safety Cabinets

- Safely store corrosives under a workbench or lab counter for quick and easy access
- These cabinets have the same quality safety features as the cabinets for flammables
- Double-wall construction, dual vents, grounding wire connections, adjustable shelves, leakproof sills, three-point self-latching doors and leveling feet
- Manual-close door(s) opens to a full 180 degrees and self-latches when pushed closed
- Superior puncture & abrasion resistance



Thomas No.	Description	Each	List	Sale
8246B12	Sure-Grip EX Undercounter Corrosives/Acid Steel Safety Cabinet, Blue, 1 shelf, 2 manual-close doors	1	1081.39	869.00

Sure-Grip® EX Corrosives/Acid Steel Safety Cabinets

Research environments using many different types of chemicals have special storage requirements. Whether using acids, bases or solvents, it is critical that incompatible liquids are segregated to avoid disastrous explosions or fire. The SDS (Safety Data Sheet) should be reviewed to determine chemical characteristics and recommended storage practices.




- Blue steel safety cabinets for corrosives have the same quality safety features as the cabinets for flammables
- Manual-close door(s) opens to a full 180 degrees and self-latches when pushed closed
- Include polyethylene trays attached to galvanized steel shelves and a separate polyethylene liner for the bottom sump resist aggressive chemicals
- Self-close door(s) shuts and latches automatically when a fusible link melts at 165°F (74°C) under fire conditions
- Manual close models meet OSHA and NFPA 30

Thomas No.	Description	Each	List	Sale
8246B16	Sure-Grip EX Undercounter Corrosives/Acid Steel Safety Cabinet, Blue, 22 gal, 1 shelf, 2 manual-close doors	1	1131.49	909.27



Polyethylene Safety Disposal Cans

Eagle Poly Disposal Safety Cans provide an intermediate station for temporary storage of liquid chemical wastes.



- Cans are FM approved and feature a wide 4" spout
- Safety cans have a cap that locks open to 80 degree angle for easy filling
- High density polyethylene (HDPE) construction resists dents and punctures
- Stainless steel fittings and flame arrestor, ideal for aggressive chemicals and corrosives

Thomas No.	Description	Each	List	Sale
1197B16	Disposal Can, Polyethylene, 3 Gallon, Red	1	261.14	167.88
1197B17	Disposal Can, Polyethylene, 3 Gallon, White	1	261.14	167.88
1197B18	Disposal Can, Polyethylene, 3 Gallon, Yellow	1	261.14	167.88
1197B19	Disposal Can, Polyethylene, 5 Gallon, Red	1	272.86	181.90
1197B20	Disposal Can, Polyethylene, 5 Gallon, White	1	282.96	181.90
1197B21	Disposal Can, Polyethylene, 5 Gallon, Yellow	1	282.96	181.90

3M

Virtua™ Protective Eyewear

- Three size options to fit a variety of facial profiles
- All stations are ANSI Z308,1-2003 compliant



Thomas No.	Description	Each	List	Sale	Case	List	Sale
1220W32	Clear Temples, Clear Hard Coat Lens	1	3.41	1.40	20	68.20	28.00
1186W35	Clear Temples, Clear Hard Coat Lens	1	4.02	2.04	—	—	—
1184L66	Clear Temples, Clear Hard Coat Lens	1	2.98	1.95	—	—	—

Honeywell

Banded Ear Plugs

- Highest attenuation among banded hearing protectors
- Soft ear pads mold quickly for a comfortable, custom fit



Thomas No.	Description	Each	List	Sale
8245A01	Banded Earplugs, NRR 2	1	9.18	4.60
8245A03	Replacement Ear Plugs, Pair	1	5.39	2.61

MICROFLEX
THE MOST TRUSTED NAME IN GLOVES

SafeGrip® Latex Gloves

- Blue, textured fingers
- Ambidextrous and Powder-free



Thomas No.	Description	Box	List	Sale	Case	List	Sale
5761W21	SafeGrip, Gloves, Small	50	24.02	15.02	500	216.20	150.20
5761W26	SafeGrip, Gloves, Medium	50	24.02	15.02	500	216.20	150.20
5761W31	SafeGrip, Gloves, Large	50	24.02	15.02	500	216.20	150.20
5761W36	SafeGrip, Gloves, X-Large	50	24.02	15.02	500	216.20	150.20

Honeywell

NORTH

First Aid Stations

- Ideal for office, lab or manufacturing settings
- All stations are ANSI Z308,1-2003 compliant



Thomas No.	Description	Each	List	Sale
1233T96	First Aid Station, 100 Person	1	184.55	131.81
1233T97	Lab Soaker, Light Duty Pad, 20" X 16'	1	265.51	189.64

Bulk First Aid Kits

- Contains economic bulk supplies that can be readily accessed when needed
- Metal kits are made of steel, plastic kits are waterproof




Thomas No.	Description	Each	List	Sale
1233T99	Bulk First Aid Kit, 10 Person, Metal	1	40.26	28.76
1233U00	Bulk First Aid Kit, 25 Person, Metal	1	42.30	30.21
1233U02	Bulk First Aid Kit, 75 Person, Metal	1	106.38	75.99
1233U03	Bulk First Aid Kit, 10 Person, Plastic	1	24.25	17.31
1233U04	Bulk First Aid Kit, 25 Person, Plastic	1	28.46	20.33
1233U05	Bulk First Aid Kit, 50 Person, Plastic	1	56.47	40.33
1233U06	Bulk First Aid Kit, 75 Person, Plastic	1	86.36	61.69

NPS CORP

Multi-Purpose Lab Soakers

- Spillyfyer® Soakers provide absorption protection against spills and cushioning to minimize breakage
- Made with multiple layers of 100% recycled tissue and a plastic poly backing for diaper-like characteristics
- Heavy duty models feature durable fabric top



Thomas No.	Description	Case	List	Sale
1233K17	Light Duty Roll, 400'L x 20"W	2	186.84	153.90
1233K18	Light Duty Pad, 20"L x 16"W	250	85.76	70.68
1233K19	Heavy Duty Roll w/ Fabric Top, 200'L x 20"W	2	151.36	124.66
1233K20	Heavy Duty Pad w/ Fabric Top, 20"L x 16"W	200	112.01	92.33

Kimberly-Clark

JACKSON SAFETY* G80 NITRILE Chemical Resistant Gloves

- Provide protection against acids, caustics, bases, solvents and other aqueous solutions
- Enhanced chemical resistance, improved durability and comfort, as well as environmentally improved packaging
- Textured fingertips and palm for enhanced grip
- Now approved for food contact in the US and Canada



Thomas No.	Description	Case	List	Sale
1202A89	JACKSON SAFETY* G80 NITRILE Chemical Resistant Gloves, Size 7 - SM	60	145.84	110.17
1202A90	JACKSON SAFETY* G80 NITRILE Chemical Resistant Gloves, Size 8 - MD	60	145.84	110.17
1202A91	JACKSON SAFETY* G80 NITRILE Chemical Resistant Gloves, Size 9 - LG	60	145.84	110.17
1202A92	JACKSON SAFETY* G80 NITRILE Chemical Resistant Gloves, Size 10 - XL	60	145.84	110.17
1202A93	JACKSON SAFETY* G80 NITRILE Chemical Resistant Gloves, Size 11 - XXL	60	145.84	110.17

slice™

Auto-Retractable Box Cutter

- Auto-retractable for extra safety
- High visibility green rubberized interior grip
- Includes one double-sided ceramic blade



Thomas No.	Description	Case	List	Sale
1191X86	Box Cutter, Auto-Retractable, Ceramic	6	203.20	139.90

Manual Pen Cutter

- Manual 3-position slider button
- Includes one double-sided ceramic blade
- Stays sharp up to 10X longer than steel



Thomas No.	Description	Case	List	Sale
1191X98	Pen Cutter, 3 Position, Manual, Ceramic	12	264.23	181.87

Order Now!

www.thomassci.com Phone 800.345.2100 Fax 800.345.5232

Use Promo Code 0217 to receive the featured pricing. Valid until 12/31/17.



Right-to-Know, Safety-Vented Wash Bottles w/GHS Labeling

- Wide mouth 53mm (2") bottle for easy refilling
- Color-coded polypropylene closures aid in quick and easy identification
- Unique valve prevents pressure build-up eliminating potentially dangerous chemical drips



Thomas No.	Description	Pack	List	Sale	Case	List	Sale
1199R27	Acetone, 500mL	4	58.00	32.78	24	348.00	196.68
1199R28	Dichloromethane, 500mL	4	52.26	46.52	—	—	—
1199R29	Deionized Water, 500mL	4	58.00	32.78	24	348.00	196.68
1199R30	Distilled Water, 500mL	4	52.26	46.52	—	—	—
1199R31	Ethyl Acetate, 500mL	4	52.26	46.52	—	—	—
1199R32	Isopropanol, 500mL	4	58.00	32.78	24	348.00	196.68
1199R33	Machine Oil, 500mL	4	52.26	46.52	—	—	—
1199R34	Methanol, 500mL	4	58.00	32.78	—	—	—
1199R35	Methyl Ethyl Ketone, Green Cap, 500mL, Non-Vented	4	52.26	46.52	—	—	—
1199R36	Isotonic Saline, 500mL	4	52.26	46.52	—	—	—
1199R37	Bleach, 500mL	4	52.26	46.52	—	—	—
1199R38	Toluene, Red Cap, 500mL, Non-Vented	4	52.26	46.52	—	—	—
1199R39	Water, 500mL	4	52.26	46.52	—	—	—
1199R40	Ethanol, 500mL	4	58.00	32.78	—	—	—
1199R41	Assortment, 500mL	4	52.26	46.52	—	—	—
1199R42	Acetone, 1000mL	4	39.00	34.71	—	—	—
1199R43	Distilled Water, 1000mL	4	39.00	34.71	—	—	—
1199R44	Isopropanol, 1000mL	4	39.00	34.71	—	—	—
1199R45	Methanol, 1000mL	4	39.00	34.71	—	—	—
1199R46	Bleach, 1000mL	4	39.00	34.71	—	—	—
1199R47	Ethanol, 1000mL	4	39.00	34.71	—	—	—
1199R48	Assortment, 1000mL	4	77.93	69.37	—	—	—

Sempermed®

SemperShield® Nitrile Gloves

- Extra thickness and an extended cuff
- Excellent chemical resistance
- Tested for use with chemotherapy drugs
- Manufactured without natural rubber
- Powder-Free



Thomas No.	Description	Box	List	Sale	Case	List	Sale
1186P92	5 mil, 12" Extended Cuff, S	50	6.74	6.49	500	67.40	64.90
1186P93	5 mil, 12" Extended Cuff, M	50	6.74	6.49	500	67.40	64.90
1186P94	5 mil, 12" Extended Cuff, L	50	6.74	6.49	500	67.40	64.90
1186P95	5 mil, 12" Extended Cuff, XL	50	6.74	6.49	500	67.40	64.90
1186P96	5 mil, 12" Extended Cuff, XXL	50	6.74	6.49	500	67.40	64.90



XCEED® Powder-Free Nitrile Exam Gloves



- **Unparalleled Comfort:** In a comprehensive study XCEED® gloves markedly reduced muscle effort
- **Improved Strength:** Up to 60% stronger than leading nitrile brands and up to 90% stronger than leading latex brands
- **Leading Protection:** Setting the standard with pinhole rates that are nearly 75% below the already stringent level set by the FDA
- **Value:** Packed 250 gloves per box (230 for X-Large), XCEED® gloves save you money, time and storage space

Traditionally, glove users have had to compromise strength in order to receive a more comfortable lighter weight glove. Less comfortable gloves can increase the muscle effort required to perform repetitive tasks. Prolonged muscle effort over longer periods of time can cause hand fatigue which has been associated with more serious injuries.

Thomas No.	Description	Box	List	Sale	Case	List	Sale
1220U55	XCEED® Powder-Free Nitrile Exam, XS	250	33.29	19.62	2500	299.61	196.20
1220U56	XCEED® Powder-Free Nitrile Exam, S	250	33.29	19.62	2500	299.61	196.20
1220U57	XCEED® Powder-Free Nitrile Exam, M	250	33.29	19.62	2500	299.61	196.20
1220U58	XCEED® Powder-Free Nitrile Exam, L	250	33.29	19.62	2500	299.61	196.20
1220U59	XCEED® Powder-Free Nitrile Exam, XL	230	33.29	19.62	2300	299.61	196.20



Simple Green Safety Towels®

- Specially designed with the industrial cleaning strength of Crystal Simple Green infused into thermally bonded, dual textured, melt-blown, lint-free polypropylene towels
- Safe and effective way to remove grease, oil, wet paint, carbon, graphite, tar, soot, grime, most inks and other soils from: metals finished wood, laminates, fiberglass, porcelain, tile, ceramics, concrete, terrazzo, carpets ... and other industrial surfaces
- Pre-moistened with non-toxic and biodegradable Crystal Simple Green®

Multi-purpose cleaning towels can handle just about any clean-up challenge.



Thomas No.	Description	Case	List	Sale
1200X47	Simple Green Safety Towels, 75 Count	6	167.39	108.94

Lemon Scent Simple Green® All-Purpose Cleaner

- Original non-toxic, biodegradable, concentrated formula is now available in lemon scent
- Can be diluted with water to remove grease, soap scum, coffee and juice stains, ink, lipstick and more from any washable, indoor or outdoor surface
- Lemon Scent is ideal for kitchens and baths!



Thomas No.	Description	Case	List	Sale
1200X70	Lemon Scent Simple Green All-Purpose Cleaner, 24 oz.	12	139.18	106.30
1200X71	Lemon Scent Simple Green All-Purpose Cleaner, 1 gal.	6	143.45	109.47

Thomas
Scientific

We Believe You Are Important
— How Can We Help?

800.345.2100 • 800.345.5232 fax
PO Box 99 • Swedesboro, NJ 08085
thomassci.com • value@thomassci.com

