



Keep a Clean Lab

with a Selection of Thomas Scientific's
Clean-Up Products!



800.345.2100 • 800.345.5232 fax
PO Box 99 • Swedesboro, NJ 08085
www.thomassci.com • value@thomassci.com



Alconox



Alconox Detergent

- For general laboratory use
- Suitable in ultrasonic cleaners
- Mild, odorless, nontoxic powdered wetting agent and detergent for cleaning glassware, porcelain, metal, plastic or rubber

Detergent removes grease, blood, tissue, etc. and is especially effective for penetrating irregular and inaccessible surfaces. Usual dilution 1 tablespoon to 1 gallon water gives pH 9.5 solution.

Thomas #	Manu #	Description
2902G05	1104	Alconox, 4-Lb. Bx.
7296B55	AG-8101	Abrasion Tester, 230 V, 50 Hz
7296B60	AG-8110	Brush Holder And Mat, 1.5 X 3.5'

Anti-Static Brushes



Thomas[®] Anti-Static Static-Away[™] Brush

- Brushes remove static charges, resulting in the adhesion of particles such as dust to an object
- Soft hair brush designs prevents scratching
- Ideal for dusting balances, cuvettes, microscope slides, optics and instrument/computer screens
- Removes static as it brushes
- Conductive composite handle
- Long handle allows access to difficult places and eliminates the possibility of a hand touching
- Brushes weigh less than an ounce

Thomas #	Description
1227T91	Brush 1" Static Away, 5" Handle
1227T92	Brush 2" Static Away, 5" Handle
1227T93	Brush 3" Static Away, 6" Handle

Cheese Cloth



Thomas® Cheesecloth Wipers

- For cleaning optics, wiping up spills or drying delicate labware
- 100% pure cotton, clean, absorbent and non-abrasive
- Will absorb six times their own weight of solvents and aqueous solutions
- Reagent grade quality
- Meets all lab requirements

Lint-free 9 inch squared 8-ply pad (2904A20) unfolds to 18 x 36 inch wipes and are packed 60 to a dispenser box. Brilliant white, 4 inch square 8-ply pad (1235C29) unfolds to approximately 7 x 14 inches and are packed 200 to a box.

Thomas #	Description
1235C29	Mini-Cheesecloth Wipes, Pack 200
2904A20	Cheesecloth Wipes, Box 60

Compressed Air



Dust-Off

- 100% Ozone-safe
- Pure moisture-free blasts
- Extension tubes for precise dusting

Ideal for cleaning in laboratories, photographic equipment and computers. The Special Application Formula (2905T03) has been formulated for use in especially sensitive environments, such as scientific research labs and other areas where flammability is a prime concern. Patented Vector Valve allows user to angle the blast 360° into places that would normally be inaccessible. Fits Dust-Off Plus models.

Thomas #	Manu #	Description
1231K54	DPVC	Dust-Off Plus Replacement Vector Valve
2905T01	DPSXL	Dust-Off XL, Disposable, 10 oz.
2905T03	DPNXL	Dust-Off XL, Disposable, Special Application Formula, 10 oz.
2905T07	DPNR	Dust-Off Plus, 8 oz. Special Application Refill For Refillable Canisters

Compressed Air



Thomas® Whoosh-Duster

- 100% Ozone-safe
- Pure moisture-free blasts
- Extension tubes for precise dusting

Ideal for cleaning in laboratories, photographic equipment and computers. The Special Application Formula (2905T03) has been formulated for use in especially sensitive environments, such as scientific research labs and other areas where flammability is a prime concern. Patented Vector Valve allows user to angle the blast 360° into places that would normally be inaccessible. Fits Dust-Off Plus models.

Thomas #	Manu #	Description
2905S44	DPVC	Whoosh-Duster, 11 oz.

Glass Disposal



Thomas® Glass Disposal Boxes

- Safely handle broken glass
- Sturdy corrugated cardboard
- 2 mil thick polypropylene liner
- Two models: bench and floor

Keep the lab safe with Thomas glass disposal boxes. Constructed from rugged corrugated cardboard and supplied with a 2 mil thick puncture-resistant polyethylene liner, the unit safely contains broken glass and spilled liquids. When full, simply pull the safety cap into place and dispose of the entire unit. Striking graphics make the glass disposal box easy to find. Box comes in two versions: Benchtop (8" x 8" x 10") and Floor (12" x 12" x 27") which hold 5.5 and 40 lbs respectively.

Thomas Glass Disposal Box Cover

- Dishwasher-safe and autoclave-safe at 112°C (250°F)
- Reusable cover opens quickly and easily by pressing the built-in release tab
- Fits any disposal carton that has a 12" x 12" opening
- Provides a wider opening (8.25 x 9.5") than cardboard lids, permitting safer and more sanitary deposit of wastes

Thomas #	Description
1236E00	Glass Disposal Box, Floor
1236E01	Glass Disposal Box, Bench
1225Z68	Glass Disposal Box Cover

Gloves

Diamond Grip Plus Latex Gloves



- Powder-free, latex, exam gloves
- Ambidextrous
- Fully textured surface
- Tapered fingertips

Non-sterile gloves are made of natural rubber latex and are fully textured for extraordinary performance. Powder-free to help reduce the risk of skin irritation or contamination from glove powders. For single use only. Cuff thickness: 5.7 ±1.0 mils; palm thickness: 7.7 ±1.0 mils; finger thickness 7.9 ±1.0 mils. Length: 9.0 inches.

Caution: This product contains natural rubber latex (latex), which may cause allergic reactions. Safe use of this glove by or on latex sensitive individuals has not been established.

Thomas #	Manu #	Description
5762E20	DGP350XS	Diamond Grip Plus Gloves, Extra Small
5762E24	DGP350S	Diamond Grip Plus Gloves, Small
5762E28	DGP350M	Diamond Grip Plus Gloves, Medium
5762E32	DGP350L	Diamond Grip Plus Gloves, Large
5762E36	DGP350XL	Diamond Grip Plus Gloves, Extra Large

Supreno SE Nitrile Gloves



- Designed for exceptional comfort, strength and durability
- No latex sensitivity concerns
- Textured fingertips
- Interior polymer coating

Supreno SE blue gloves are designed to provide the chemical compatibility of nitrile without fatigue during prolonged use. Textured fingers allow for improved grip in both wet and dry applications. Powder-free gloves avoid the mess inherent of powdered gloves by using a polymer coating to allow for easier donning and doffing. Supreno SE gloves meet or exceed ASTM specifications for thickness, tensile strength and elasticity. Ambidextrous gloves are 9.4 inches long with a beaded cuff. Thickness is 5 mils at the cuff, 6.3 mils in the palm and 7.1 mils in the textured fingertips.

Thomas #	Manu #	Description
5761Q31	SU690XS	Supreno SE Gloves, Extra Small
5761Q35	SU690S	Supreno SE Gloves, Small
5761Q39	SU690M	Supreno SE Gloves, Medium
5761Q43	SU690L	Supreno SE Gloves, Large
5761Q47	SU690XL	Supreno SE Gloves, Extra Large

Spill Kits

Spillsolv Chemical Spill Kits



- Multi-Purpose Spill Kit
- High capacity neutralizers
- Wall-mountable or portable case
- Shaker-top containers
- Neutralization indicated by color change

Uses less neutralizer which provides faster clean-ups and less expense. Forms into small clumps making clean-up and disposal quicker and simpler. Sturdy, bright yellow case is lightweight and resists deterioration. Shaker-top bottles eliminate having to scoop neutralizers out of a container. Replacement kits contain neutralizers only.

- For acid, caustic and solvent spills
- For multi-purpose, acid, solvent and caustic kits

Thomas #	Manu #	Grade	Pkg. Size	Pkg. Type
8237A01	SX1300-1	SpillSolv	1PC	Poly Box
8237A03	SX1300R-1	SpillSolv	1PAC	Carton
8237A05	SX1310-1	SpillSolv	1PAC	Poly Box
8237A07	SX1310R-1	SpillSolv	1PAC	Carton
8237A09	SX1330-1	SpillSolv	1PC	Poly Box
8237A11	SX1330R-1	SpillSolv	1PC	Carton
8237A13	SX1320-1	SpillSolv	1PC	Poly Box
8237A15	SX1320R-1	SpillSolv	1PC	Carton
8237A17	SX1300S-1	SpillSolv	1PAC	Carton
8237A19	SX1340-1	SpillSolv	1PC	Poly Box
8237A21	SX1340R-1	SpillSolv	1PC	Carton
8237A23	SX1340S-1	SpillSolv	1PC	Carton
8237A25	SX1350-1	SpillSolv	1PC	Poly Box

Wipers



Thomas® Clean-Wipes

- 70% isopropyl or 100% pure water treated wipes
- High-quality, ultra-clean grade
- Strong, soft and lint-free

Pre-moistened Clean-Wipes are just right for cleaning fragile equipment and clean up even the dirtiest parts. Non-woven cloths are moistened with either a 70/30 blend of pure, reagent-grade isopropyl alcohol and reagent-grade deionized water or 100% pure reagent-grade deionized water. Dry wipes are available. User can add any solution. Perfect for cleaning optical parts, fragile glassware, electrodes, cuvettes, microscopes or lenses.



Instantly clean metal, plastic, glass, rubber and epoxies. May be used for cleaning electronics and probes. Designed for all applications in environmentally sensitive areas. Very strong, cleanroom-grade, nonwoven cloth is soft and lint-free. Constructed of cellulose and polyester, they will not scratch delicate surfaces. They meet all lab and production requirements for cleanliness, softness, absorbency and strength. Long-lasting wipes never fall apart like paper towels.

Wipes are individually dispensed and ready for instant use, with a unique opening that allows wipes to virtually glide out when pulled. Reclosable benchtop container assures a clean wipe each time and the snap-on sealed top eliminates contamination from dust or dirt. Each canister is shrink-wrapped to maintain cleanliness during shipping. Clean-Wipes are 6 x 9 inches and packaged 100 per polyethylene canister.

Thomas #	Description
2903J75	Clean-Wipes, Deionized Water, Can 100
2903J78	Clean-Wipes, 70% Iso Alcohol, Can 100
2903J81	Clean Wipes, Dry, Can 100



Versi-Dry Light Duty Roll, 20" X 300'

Versi-Dry Soakers speed clean-up time and protect lab table tops from acids, dyes, and radioactive applications. Soakers feature THIRSTY CELLS* that quickly absorb spills. Soft surface cushions breakable equipment. Nonskid, chemical-resistant, and waterproof polyethylene backing holds soaker securely in place.

Thomas #	Manu #	Description
1228B07	3706600	Versi-Dry Light Duty Roll, 20" X 300'

Wipers



Versi-Dry Soakers

Versi-Dry Soakers speed clean-up time and protect lab table tops from acids, dyes, and radioactive applications. Soakers feature THIRSTY CELLS* that quickly absorb spills. Soft surface cushions breakable equipment. Nonskid, chemical-resistant, and waterproof polyethylene backing holds soaker securely in place.

Thomas #	Manu #	Description
1228B09	7420000	Versi-Dry Super Absorbant Floor Roll 40" X 200'
1228B10	7421800	Versi-Dry Super Absorbant Fume Hood Mat 18" X 20"
1228B11	7424000	Versi-Dry Super Absorbant Floor Mat 18" X 40"