

FREE
Lab Equipment
with purchase of item below



Value: \$279.00

8 x 1.5/2.0ml
or 4 x PCR Strips



Value: \$200.84

2,800rpm (fixed speed)



Value:
\$345.04

Accepts tubes 1.5 to 50ml

Your choice of: **MiniMouse™ II** Item: 1154J27 **Mini Vortexer** Item: 1154J94 **BioRocker™ FX** Item: 1154J74



Up to 70°C

INSTANT Flask Attachment
Patented

labForce Mini Shaking Incubator

- Touch screen control of time, temperature and shaking speed
- **MAGic Clamp™** accessories for instant magnetic clamp exchange (patented)
- **NEW!** Adjustable clamp for 125, 250 and 500ml flasks
- **Flask Capacity:** 12 x 125ml, 8 x 250ml, 5 x 500ml or 2 x 1L

NEW!
Adjustable Clamp



| | | |
|---------|---|-----------|
| 1141U52 | Mini Shaking Incubator with FREE MiniMouse Centrifuge | \$2934.59 |
| 1141U53 | Mini Shaking Incubator with FREE Mini Vortexer | \$2934.59 |
| 1141U54 | Mini Shaking Incubator with FREE BioRocker FX | \$2934.59 |
| 1222U14 | MAGic Clamp universal platform for flasks & tube racks | \$254.58 |
| 1177X57 | MAGic Clamp Adjustable, for 125ml, 250ml & 500ml flasks | \$41.41 |
| 1222U28 | MAGic Clamp, for 1L flask (max. 2) | \$60.48 |
| 1222U27 | MAGic Clamp, for 1 x microplate (max. 4) | \$79.09 |
| 1222U18 | MAGic Clamp tube rack, 30x15ml & 20x50ml tubes (max. 1) | \$157.38 |



labForce BioRocker™ 2D Rocker

- Adjustable speed & tilt for a wide range of applications
- Large 14x 12 in. work surface
- Cold room and incubator safe

| | | |
|---------|--|----------|
| 1141U55 | BioRocker 2D Rocker with FREE MiniMouse Centrifuge | \$771.44 |
| 1141U56 | BioRocker 2D Rocker with FREE Mini Vortexer | \$771.44 |
| 1141U57 | BioRocker 2D Rocker with FREE BioRocker FX | \$771.44 |

labForce High Speed Microcentrifuges



Specifications:

| | |
|----------------|----------------------------|
| Speed Range: | 500rpm~15,000rpm |
| Max. RFC: | 21,4000 xg |
| Max. Capacity: | 24x1.5/2.0ml |
| Temp. Range: | -10 to 40°C (1134Q52 only) |
| Timer: | 1s-99 min or Short-spin |

- High Speed, 21,400xg
- Quiet and compact
- Advanced, full color display
- Available with or without refrigeration

| | | |
|---------|--|-----------|
| 1141U58 | 24-M Microcentrifuge with FREE MiniMouse Centrifuge | \$2843.39 |
| 1141U59 | 24-M Microcentrifuge with FREE Mini Vortexer | \$2843.39 |
| 1141U60 | 24-M Microcentrifuge with FREE BioRocker FX | \$2843.39 |
| 1141U61 | 24-MR Microcentrifuge with FREE MiniMouse Centrifuge | \$6590.44 |
| 1141U62 | 24-MR Microcentrifuge with FREE Mini Vortexer | \$6590.44 |
| 1141U63 | 24-MR Microcentrifuge with FREE BioRocker FX | \$6590.44 |

labForce PROMO



MyTemp™ Mini Digital Incubator

- Personal sized, 13 x 14.5in. footprint (20L)
- Accepts bottles and flasks up to 2L
- Heats & Cools: Ambient -15°C to 60°C

Ordering Information:

Thomas No. 1154J61
Mini Digital Incubator \$746.78

Thomas No. 1144T88
Buy 1 Get 1 Price \$845.78 (2 ea.)

Buy 1, Get 1 for
\$99

Digital Hotplate Stirrer

- Large, Digital display
- Safety LED indicates hot surface
- Control actual temperature with optional probe

Max Temp: 380°C Max Speed: 1500rpm

Ordering Information:

Thomas No. 1154K00
Digital Hotplate Stirrer \$380.09

Thomas No. 1144T91
Buy 1 Get 1 Price \$479.09 (2 ea.)

Buy 1... Get 1 for
\$99



IncuBlock™ Plus Digital DryBaths

- Hinged cover improves uniformity and conserves energy
- Optional external temperature probe for control of actual sample temperature
- Unique Quick-Flip™ universal micro-block

Ordering Information:

Thomas No. 1154J44
labForce IncuBlock™ Plus, Single \$410.97

Thomas No. 1154J46
labForce IncuBlock™ Plus, Dual \$513.97

Thomas No. 1227U85
Block 24x1.5/2.0ml Centrifuge tubes \$137.59

Thomas No. 1227U89
Block 20x13mm Test tubes \$137.59

Thomas No. 1227U87
Block 12x15ml Centrifuge tubes \$137.59

Thomas No. 1227U90
Block 5x50ml Centrifuge tubes \$137.59

Thomas No. 1203F36
Quick-Flip™ Block, 24x1.5ml tubes, or 32x0.2ml and 14x0.5ml tubes \$171.26



Specifications:

Temp. Range: Ambient +5° to 105°C
Temp. Accuracy: +0.1°C
Temp. Uniformity: +/- 0.3°C
Temp. Increments: 0.1°C
Timer: 1 min. to 48 hrs. / continuous



Vortexer V-32 Vortex Mixer

- Unique Multi-Head™ for general purpose vortexing PLUS holds up to eight microtubes
- Powerful motor for INSTANT vortexing
- Continuous or "touch" operation
- Q-Drive™ dynamic balancing system
- Variable speed form 200 to 3400 rpm

Specifications:

Speed Range: 200 to 3400rpm
Operating Modes: Touch or continuous
Orbit: 3mm
Temp. Range: +4 to 65°C

Ordering Information:

Thomas No. 1154J84
labForce Vortexer V-32 \$283.24



MicroPlate Centrifuge

- Quickly spins down droplets in PCR plates
- Less than half the size of traditional centrifuges
- Unique, patented rotor design prevents sample spillage
- Accepts all popular PCR plates (skirted, semi and non) and most microtiter plates

Specifications:

Speed Range: 2,550 rpm / 600 xg
Capacity: 2 x 2PCR Plates
2 x microtiter plates (up to 25mm)
24 x 0.2ml PCR strips
Deceleration: 4 seconds

Ordering Information:

Thomas No. 1154J32
labForce Microplate Centrifuge \$689.00

