



50% OFF

Fischer Technical Company

E1000-05 VARIABLE VOLUME PIPETTOR

Fischer Technical variable volume pipettors are adjustable between 500-5000 μ L and feature an analog volume adjustment. To adjust the volume, simply turn the plunger until the graduations line up with your desired volume. Accuracy is $\pm 1.75\%$ and repeatability is $\pm 0.4\%$. Accepts standard tips including the Fischer Technical D893 tips (Thomas No: 1200D25).

Shipping weight 1.4lbs.

Thomas No.	Manuf. No.	Description	Quantity	List Price	Promo Price*
1200D22	E1000-05	500-5000ul adjustable micropipettor	Each	\$96.46	\$48.23

*While supplies last

Disclaimer: Pricing and offers subject to change without notice.