



See **Special Buy One, Get One**
Offers on **pages 4 and 8**

Consumables *and* Equipment Brochure

labForce
powered by *Thomas Scientific*

labForce Specials

Use promo code **0120** to
receive your special pricing.

Prices good through April 30, 2020

Ultra EZpette™

Single & Multi Channel High Precision Digital Pipettors

Ultra EZpette™ provides both accurate and precise liquid dispensing in a comfortable ergonomic design. Fully autoclavable, UV resistant.

Features

- Dual volume setting system
- Volume lock
- Easy calibration
- Adjustable ejector — **NEW “click-to-fix” system!**
- Durable sealing system effective for 60,000 cycles



Thomas No.	Description	UOM	List Price	Sale Price
1159Z35	Single Channel Blue, 0.2-2µL	1/Each	\$392.40	\$157.00
1159Z29	Single Channel Blue, 0.5-10µL	1/Each	\$392.40	\$157.00
1159Z33	Single Channel Blue, 2-20µL	1/Each	\$392.40	\$157.00
1159Z27	Single Channel Blue, 10-100µL	1/Each	\$392.40	\$157.00
1159Z31	Single Channel Blue, 20-200µL	1/Each	\$392.40	\$157.00
1159Z25	Single Channel Blue, 100-1000µL	1/Each	\$392.40	\$157.00
1159Z37	Single Channel Blue, 1000-5000µL	1/Each	\$392.40	\$157.00
1159Z23	Single Channel Blue, 1000-10000µL	1/Each	\$392.40	\$157.00

Thomas No.	Description	UOM	List Price	Sale Price
1159Z40	8-Channel, 1-10µL	1/Each	\$676.88	\$372.30
1159Z43	8-Channel, 5-50µL	1/Each	\$676.88	\$372.30
1159Z41	8-Channel, 20-200µL	1/Each	\$676.88	\$372.30
1159Z42	8-Channel, 50-300µL	1/Each	\$676.88	\$372.30
1159Z44	12-Channel, 1-10µL	1/Each	\$780.55	\$390.30
1159Z47	12-Channel, 5-50µL	1/Each	\$780.55	\$390.30
1159Z45	12-Channel, 20-200µL	1/Each	\$780.55	\$390.30
1159Z46	12-Channel, 50-300µL	1/Each	\$780.55	\$390.30

labForce Sharp Elite™ Pipette Tips

Tip features include:

- Ultralow retention resin maximizes your sample recovery
- Ultraclear resin allows sample volume to easily be visualized through the tips
- Graduated for quick confirmation of sample volumes
- 10µL XL tip features graduations at: 1, 2 and 10µL for accurate low-volume pipetting
- 1000µL XL tip features extra headspace allowing for full 1mL draws without risk of contacting the filter
- Certified free of detectable RNase, DNase, DNA, PCR inhibitors, and pyrogen-free
- Multi-channel pipettor friendly




Sharp Elite Pipet Tips provide the next generation of precision and convenience. Designed to precisely fit most universal pipettors, Sharp Elite tips deliver unmatched accuracy for your research needs. With a full range of options (bulk, racked, racked sterile, barrier/filter tips, and reloads).

Barrier Tips, Sterile, Racked

Thomas No.	Description	UOM	List Price	Sale Price
1149J69	10µL Barrier Tips, Racked	960/Pack	\$78.36	\$50.95
1149J70	10µL XL Barrier Tips, Racked	960/Pack	\$78.36	\$50.95
1149J71	20µL Barrier Tips, Racked	960/Pack	\$78.36	\$50.95
1149J88	100µL Barrier Tips, Racked	960/Pack	\$78.36	\$50.95
1149J72	200µL Barrier Tips, Racked	960/Pack	\$78.36	\$50.95
1138W95	200µL XL Barrier Tips, Racked	192/Pack	\$23.78	\$15.50
1149J73	300µL Barrier Tips, Racked	960/Pack	\$85.36	\$55.50
1145N16	1000µL XL Barrier Tips, Racked	960/Pack	\$84.92	\$55.20

Color coded wafers and big bold front of rack labels allow for quick tip volume identification.



Racked Tips

Thomas No.	Description	UOM	List Price	Sale Price
1149J90	10µL Racked Tips	960/Pack	\$44.45	\$28.90
1149J93	10µL XL Racked Tips	960/Pack	\$45.95	\$29.90
1149K10	200µL Racked Tips	960/Pack	\$42.94	\$27.90
1138M55	200µL XL Racked Tips	192/Pack	\$12.32	\$8.05
1138M56	300µL Racked Tips	960/Pack	\$59.44	\$38.65
1145N14	1000µL Racked Tips	960/Pack	\$55.90	\$36.35

Bulk Tips

Thomas No.	Description	UOM	List Price	Sale Price
1149J97	10µL Bulk Tips	1000/Pack	\$22.52	\$14.65
1148B39	10µL XL Bulk Tips	1000/Pack	\$21.98	\$14.30
1149J98	200µL Bulk Tips	1000/Pack	\$15.95	\$10.40
1148F60	300µL Bulk Tips	1000/Pack	\$30.84	\$20.05
1145N12	1000µL XL Bulk Tips	1000/Pack	\$32.25	\$21.00



Sterile, Racked Tips

Thomas No.	Description	UOM	List Price	Sale Price
1149J91	10µL Racked Tips	960/Pack	\$50.77	\$33.05
1148B42	10µL XL Racked Tips	960/Pack	\$50.77	\$33.05
1148F34	200µL Racked Tips	960/Pack	\$48.93	\$31.85
1138W98	200µL XL Racked Tips	192/Pack	\$13.65	\$8.90
1148B83	300µL Racked Tips	960/Pack	\$65.18	\$42.40
1145N15	1000µL Racked Tips	960/Pack	\$62.60	\$40.70

Bench Pack Reloads

Thomas No.	Description	UOM	List Price	Sale Price
1138M50	10µL Tips	960/Pack	\$40.52	\$26.35
1138M51	10µL XL Tips	960/Pack	\$35.70	\$23.25
1138M52	200µL Tips	960/Pack	\$33.07	\$21.50
1138M53	300µL Tips	960/Pack	\$52.28	\$34.00
1145N13	1000µL XL Tips	768/Pack	\$38.02	\$24.75

All Sharp Elite Reloads are compatible with the Sharp Elite Empty Racks.



labForce 6-Place Pipette Stand



labForce six place pipette holder that accepts most major brands of pipettes. The pipette stand is made of durable, 3/16" acrylic for years of use. The skid resistant feet will prevent the stand from sliding on your bench top.
Dimension: 11"W x 9.25"H x 6"D

Thomas No.	Description	UOM	List Price	Sale Price
1140H35	Pipette Stand, 6-Place, Clear Acrylic	1/Each	\$65.00	\$42.25

Denville™ FASTPETTE™ Pipette Controller

labForce Controller — Lightweight and ergonomically shaped motorized pipette controller designed for work with glass or plastic serological pipettes from 0.5 to 100mL range. Conveniently positioned switches allow choosing different operation modes and speed depending on the volume of the serological pipette and liquid viscosity. Additionally aspirating and dispensing speed is controlled through two comfortable finger triggers.



LCD Display Pipette

When fully charged, the NiMH batteries enable up to eight hours of continuous use. The large LCD display the serological high viscosity faspette indicates the battery status along with the pipetting mode and speed selected. In order to protect the device against overfilling, the controller is equipped with a 0.2µm PTFE filter. The filter, the pipette holder and the nosepiece are autoclavable. Also supplied with a universal battery charger, an additional hydrophobic 0.2 µm filter and a two-position charging stand.


Thomas No.	Description	UOM	List Price	Sale Price
1159Z16	labForce Pipette Controller	1/Each	\$315.90	\$221.15

Use promo code
0120
to receive your
special pricing.

labForce Serological Pipettes

labForce disposable serological pipettes with cotton plugs are ideal for tissue culture applications. Crisp, black, highly legible permanent graduations prevent ambiguity and guarantee precise reading and accurate dispensing within $\pm 2\%$. Reverse graduations show volume removed or remaining on all pipettes larger than 1mL. Generous negative graduations provide a additional pipetting volume. labForce serological pipettes are made from transparent USP Class VI polystyrene sterilized by gamma irradiation. Each lot is free from pyrogens (≤ 0.125 EU/mL), DNases, and RNases.

- Color-coded for easy Identification
- Easy-to-open paper/plastic package
- Generous negative graduation
- Individually wrapped and sterile
- Accuracy within $\pm 2\%$
- RNase/DNase, pyrogen free



Thomas No.	Description	UOM	List Price	Sale Price
1154Q54	labForce Serological Pipettes, 1mL, Sterile, Cotton Plug	100/Pack	\$49.09	\$17.20
1154Q55	labForce Serological Pipettes, 2mL, Sterile, Cotton Plug	100/Pack	\$54.84	\$19.20
1158R02	labForce Serological Pipettes, 5mL, Sterile, Cotton Plug	200/Case	\$128.05	\$44.85
1158R03	labForce Serological Pipettes, 10mL, Sterile, Cotton Plug	200/Case	\$136.44	\$47.80
1158R04	labForce Serological Pipettes, 25mL, Sterile, Cotton Plug	200/Case	\$284.25	\$99.50
1154Q56	labForce Serological Pipettes, 50mL, Sterile, Cotton Plug	100/Case	\$432.96	\$151.55

labForce MiniMouse™ II Mini Centrifuge



- *Clik-n-Lok™* Instant rotor attachment
- 8 position microtube rotor & 4 position PCR strip rotor
- Five "hot" colors
- Storage compartment for PCR rotor & 0.5mL tubes
- Includes round & PCR rotors plus tube adapters
- Near silent operation
- Decelerates in just 1 second
- Conserves valuable bench space

Technical Data

- **Speed:** 6,000rpm/2,000xg
- **Capacity:** 8 x 1.5/2.0mL tubes 32 x 0.2mL PCR tubes
- **Dimensions:** 11.4 x 1.5 x 11.4cm (4.5 x 5.9 x 4.5")
- **Weight:** 1.4kg (3 lbs) Universal Power: 100-240V, 50-60Hz

Thomas No.	Description	UOM	Price
1154J27	labForce Mini Centrifuge, Red, 100-240V (US)	1/Each	\$289.00

labForce 24-M & 24-MR High Speed Microcentrifuges

Receive a FREE labForce Digital Hot Plate Stirrer when you purchase a 115V labForce Microcentrifuge. Offer valid through April 30, 2020.



- High-speed 21,400 xg
- Quiet and compact, 11" wide
- Fast acceleration & deceleration, 14 seconds
- Additional optional rotors, 12 x 5mL and 4 x 8-PCR Strip
- LCD full color display
- 2-Year Manufacturer's Warranty
- New, convenient Snap-On rotor lid

Technical Data

- **Speed:** 15,000 rpm/21,400 xg
- **Capacity:** 24 x 1.5/2.0mL tubes; 12 x 5mL tubes
- **Temp. Range:** -10°C to 40°C (24 MR only)

FREE labForce Digital Hot Plate Stirrer!

Thomas No.	Description	UOM	Price
20A00C737	PROMOTION! 1154Q50 labForce Microcentrifuge with FREE 1140K86 Digital Hotplate Stirrer	1/Each	\$2,646.96
20A00C738	PROMOTION! 1154Q52 labForce refrigerated Microcentrifuge with FREE 1140K86 Digital Hotplate Stirrer	1/Each	\$6,364.35

labForce Amplifyt™ PCR 8-Strip Tube & Strip Cap Packs

- Thin Wall PCR Strip Tubes & Caps - 8 Strip
- Available with dome caps or optically clear flat caps for real-time PCR
- Available 0.2mL dome or 0.1mL low profile / "Fast" thermal cyclers (1149Z81)
- Available in clear, white, or assorted
- Convenient combo pack of tubes & caps – no need to purchase separately
- Deep seating caps provide a reliable, consistent seal to eliminate evaporation
- Thin walls allow for efficient heat transfer and superior yields
- Certified RNase, DNase & Pyrogen-free



Thomas No.	Description	UOM	List Price	Sale Price
1149Z65	0.2mL 8-Strip PCR Tubes & Domed Strip Caps Packed Together, Clear	125/Pack	\$95.73	\$57.45
1139Y43	0.2mL 8-Strip PCR Tubes & Domed Strip Caps Packed Together, Assorted Colors	125/Pack	\$99.56	\$59.75
1149Z81	0.1mL 8-Strip Low-Profile PCR Tubes & Flat Strip Caps Packed Together, Clear	125/Pack	\$95.73	\$57.45
1149K07	0.2mL 8-Strip PCR Tubes & Flat Strip Caps Packed Together, Clear	125/Pack	\$95.73	\$57.45
1149Z78	0.2mL 8-Strip PCR Tubes & Flat Strip Caps Packed Together, White	125/Pack	\$100.84	\$60.55

labForce Screw Cap Centrifuge Tubes



- Printed Graduations (with marking area)
- New clearer polypropylene copolymer
- High centrifugal strength - up to 15,000 RCF
- Flat markable cap with leakproof seal
- Phenol / Chloroform resistant
- DNase, RNase and pyrogen free

Thomas No.	Description	Case	List Price	Sale Price
1158R11	15mL Tubes w/Screw Cap, Sterile, Bag/25	500	\$109.68	\$82.30
1159M36	15mL Tubes w/Screw Cap, Sterile Tray/50	500	\$110.91	\$83.20
1158R09	50mL Tubes w/Screw Cap, Sterile Bag/25	500	\$138.33	\$103.75
1158R10	50mL Tubes w/Screw Cap, Sterile Tray/25	500	\$143.79	\$107.85

Use promo code
0120
to receive your
special pricing.

labForce Amplifyt™ PCR Split-Strip Tubes with Attached Caps

8 Strip Split-Strip with attached individual caps. Available with both dome caps and optically-clear flat caps for real-time. Available in two sizes to fit virtually every thermal cycler or real-time instrument. 0.2mL for standard blocks and 0.1mL for low-profile or "fast" blocks.



- Available in clear or assorted
- All caps are flat and optically-clear for realtime PCR assays
- Split-a-part design is easy to break allowing you to use only the tubes that you need
- Conveniently stacked in bags of 12 strips per bag (10 bags per unit) minimizes the risk of contamination
- Certified free of detectable RNase DNase, DNA, PCR inhibitors, and tested pyrogen-free

Thomas No.	Description	Case	List Price	Sale Price
1148A16	0.2mL Split-strip 8-Strip tubes w/individually attached dome cap, Clear	120/Pack	\$79.83	\$47.90
1148A24	0.2mL Split-strip 8-Strip tubes w/individually attached dome cap, Assorted	120/Pack	\$86.76	\$52.10
1148A25	0.2mL Split-strip 8-Strip tubes w/individually attached flat caps, Clear	120/Pack	\$84.74	\$50.85
1148A26	0.1mL Split-strip 8-Strip tubes w/individually attached flat caps, Low-profile, Clear	120/Pack	\$94.56	\$56.75

Amplifyt™

Rigid-Strip PCR Tubes with Attached Caps

8 Strip Rigid Strips with attached individual caps. Available in two sizes to fit virtually every thermal cycler or real-time instrument. 0.2mL for standard blocks and 0.1mL for low-profile or "fast" blocks.



- Available in clear, assorted, or white
- All caps are flat and optically-clear for realtime PCR assays
- Split-a-part design is easy to break allowing you to use only the tubes that you need
- Rigid design is easy to work with minimizing the need to work with a PCR rack
- Conveniently stacked in bags of 12 strips per bag (10 bags per unit) minimizes the risk of contamination
- Certified free of detectable RNase DNase, DNA, PCR inhibitors, and tested pyrogen-free

Thomas No.	Description	UOM	List Price	Sale Price
1148A28	0.2mL Rigid Strip 8-Strip tubes optically clear flat caps, Clear	120/Pack	\$80.14	\$48.10
1148A34	0.2mL Rigid Strip 8-Strip tubes optically clear flat caps, Assorted*	120/Pack	\$86.76	\$52.10
1148A35	0.2mL Rigid Strip 8-Strip tubes optically clear flat caps, White	120/Pack	\$87.77	\$52.70
1159V39	0.1mL Rigid Strip 8-Strip tubes optically clear flat caps, Clear	120/Pack	\$87.18	\$52.35
1148A36	0.1mL Rigid Strip 8-Strip tubes optically clear flat caps, White	120/Pack	\$91.23	\$54.75

*Blue, Green, Red, Violet, Yellow

Amplifyt™

96-Well Semi-Skirted PCR Plates

Broad thermal cycler compatibility, high-contrast, black alphanumeric indexing for easy well identification.

- Crystal-clear wells to easily verify sample volume samples easy
- Ultra-thin, uniform wells ensure optimal heat transfer for high reaction efficiency
- Well ridges reduce sample evaporation when sealing with film or foil
- Manufactured from high-clarity, ultra-pure polypropylene (meets FDA 21 CFR 177.1520 & USP Class VI)
- Certified free of detectable DNase, RNase, DNA, PCR inhibitors, and tested pyrogen-free



Thomas No.	Description	UOM	List Price	Sale Price
1148A75	96-Well PCR Plates, Semi-Skirt, Standard, Clear	10/Pack	\$31.72	\$19.05
1149K03	96-Well PCR Plates, Semi-Skirt, Standard, Clear	25/Pack	\$79.29	\$47.60
1148A73	96-Well PCR Plates, Semi-Skirted, Raised Rim, FAST, Low-Profile, Clear	10/Pack	\$31.72	\$19.05
1139Y45	96-Well PCR Plates, Semi-Skirted, ABI-Type, Clear	25/Pack	\$79.29	\$47.60
1148A58	96-Well PCR Plates, Semi-Skirted, ABI-Type, Clear	10/Pack	\$31.72	\$19.05

labForce

MicroPlate Centrifuges



- Quickly spins down droplets in PCR plates
- Unique rotor design prevents sample spillage (patent pending)
- 50% smaller than traditional plate centrifuges
- Accepts all popular PCR plates (skirted, semi & non) & most microtiter plates

Technical Data:

- Speed: 100rpm/400xg
- Capacity: 2 x PCR Plates; 2 x microtiter plates; 24 x 0.2mL PCR strips
- Dimensions: 23 x 26 x 19.7cm (9.2 x 10.2 x 7.75")
- Weight: 4.1kg (9 lbs)
- Universal Power: 120-230V, 50-60Hz

Thomas No.	Description	UOM	Price
1154J32	labForce MicroPlate MicroCentrifuge, 115V	1/Each	\$715.00


Use promo code
0120
to receive your
special pricing.

Amplifyt™

96-Well PCR Plates - No-Skirt

Broad thermal cycler compatibility, high-contrast, black alphanumeric indexing for easy well identification.

- 96 x 0.2mL wells or 96 x 0.1mL wells
- Flat top for ideal sealing w/automated & manual seals
- Slightly raised rim around each well for improved sealing with strip caps, films, & foils
- Cuttable version (1148A46) is easy to cut to the number of wells you need
- Certified DNase, RNase, and Pyrogen Free



Thomas No.	Description	UOM	List Price	Sale Price
1149K04	Flat Top, Standard, Clear	25/Pack	\$79.29	\$47.60
1148A38	Flat Top, Standard, Clear	10/Pack	\$31.72	\$19.05
1148A46	Non-Skirted, Cuttable, Clear	10/Pack	\$36.72	\$22.05
1148A47	Low-Profile, Non-Skirted, Cuttable, Clear	20/Pack	\$66.46	\$39.90

Amplifyt™ 96-Well PCR Full Skirt Plates

- 96-Well PCR plate, full skirt, clear or white
- 96 x 0.1mL wells (low-profile)
- Compatible with a broad range of thermal cyclers including many Biorad and Eppendorf instruments
- Flat Top for ideal sealing with automated and manual seals
- Slightly raised rim around each well for improved sealing with strip caps, films, and foils
- Vivid black alphanumerics
- Single notch design (A12)



Thomas No.	Description	UOM	List Price	Sale Price
1149K05	96-Well PCR Plates, Full-Skirted, Clear	25/Pack	\$83.07	\$49.85
1148A98	96-Well PCR Plates, Full-Skirted, Clear	10/Pack	\$33.23	\$19.95
1148B05	96-Well PCR Plates, Full-Skirted, White	10/Pack	\$39.20	\$23.55
1148B06	96-Well PCR Plates, Full-Skirted, Clear, Barcoded	10/Pack	\$39.92	\$24.00

labForce Standard Weighing Boats

- Manufactured and packaged in a clean room
- Select white or blue for best color contrast to ensure sample is completely removed
- Biologically inert and resistant to dilute and weak acids, aqueous solutions, alcohols, and bases
- Temperature range -10 to 70°C (14 to 158°F)
- Boats are shallow with a wide, flat bottom to resist tipping
- Rounded corners simplify transfer
- Anti-static versions available
- 2 colors, and 3 sizes available
- All are FDA/food-grade compliant



Thomas No.	Description	UOM	List Price	Sale Price
1163Z53	Weighing Boat, Regular, White, Small (46 x 46 x 8mm)	500/Pack	\$38.93	\$23.40
1163Z54	Weighing Boat, Anti-Static, White, Small (46 x 46 x 8mm)	500/Pack	\$41.37	\$24.85
1163Z55	Weighing Boat, Regular, White, Medium (85 x 85 x 24mm)	500/Pack	\$51.07	\$30.65
1163Z56	Weighing Boat, Anti-Static, White, Medium (85 x 85 x 24mm)	500/Pack	\$54.58	\$32.75
1163Z57	Weighing Boat, Regular, White, Large (140 x 140 x 22mm)	500/Pack	\$104.60	\$62.80
1163Z58	Weighing Boat, Anti-Static, White, Large (140 x 140 x 22mm)	500/Pack	\$110.41	\$66.25
1163Z67	Weighing Boat, Regular, Blue, Small (46 x 46 x 8mm)	500/Pack	\$39.67	\$23.85
1163Z68	Weighing Boat, Regular, Blue, Medium (85 x 85 x 24mm)	500/Pack	\$52.09	\$31.30
1163Z69	Weighing Boat, Regular, Blue, Large (140 x 140 x 22mm)	500/Pack	\$106.70	\$64.05

Thomas® Reagent Reservoirs

Choose from four styles - 25mL Divided, 25mL, 50mL and 100mL. 25mL divided reservoirs are divided into 3 separate 5mL channels; great for low volume pipetting. Sterile and certified versions are certified RNase and DNase-free and non-pyrogenic. Bulk packed reservoirs come in packs of 5 reservoirs.



- "V" bottom design allows for maximum recovery of reagents
- Graduated with convenient pouring spouts at all four corners
- Extra wide base for increased stability
- Unique "Nesting" feature makes separation of bulk packed reservoirs easy

Thomas No.	Description	UOM	List Price	Sale Price
1228K15	Reagent Reservoir, 25mL Divided, Sterile, Individually Wrapped	100/Case	\$165.25	\$82.65
1228K16	Reagent Reservoir, 25mL Divided, Sterile, Bulk Packed	200/Case	\$163.39	\$106.25
1228K17	Reagent Reservoir, 25mL Divided, Non-Sterile, Bulk Packed	100/Case	\$78.36	\$50.95
1228K18	Reagent Reservoir, 25mL, Sterile, Individually Wrapped	100/Case	\$177.14	\$88.60
1228K19	Reagent Reservoir, 25mL, Sterile, Bulk Packed	200/Case	\$141.19	\$91.80
1228K20	Reagent Reservoir, 25mL, Non-Sterile, Bulk Packed	100/Case	\$80.84	\$52.55
1228K21	Reagent Reservoir, 50mL, Sterile, Individually Wrapped	100/Case	\$207.00	\$103.50
1228K22	Reagent Reservoir, 50mL, Sterile, Bulk Packed	200/Case	\$144.88	\$94.20
1228K23	Reagent Reservoir, 50mL, Non-Sterile, Bulk Packed	100/Case	\$71.97	\$46.80
1228K24	Reagent Reservoir, 100mL, Sterile, Individually Wrapped	100/Case	\$228.22	\$114.15
1228K25	Reagent Reservoir, 100mL, Sterile, Bulk Packed	200/Case	\$156.70	\$101.90
1228K26	Reagent Reservoir, 100mL, Non-Sterile, Bulk Packed	100/Case	\$122.03	\$61.05


Use promo code
0120
to receive your
special pricing.

labForce Analytical Laboratory Balances

Standard weighing, highly precise measurement

A 3-door draft shield is standard on all analytical models (excludes milligram and high-capacity models). All draft shield doors, stainless steel weigh pan, and draft chamber floor are easy to remove for cleaning.

- Precise measurement & repeatability
- Easy-to-read information arrangement
- USB & RS-232 Interface
- User-friendly draft shield
- Five year manufacturer's warranty included on all models



Thomas No.	Description	UOM	List Price	Sale Price
1160T96	Balance, 62g x 0.1mg, 110-240V	1/Each	\$1,990.97	\$1,493.25
1160T97	Balance, 160g x 0.1mg, 110-240V	1/Each	\$2,212.80	\$1,659.60
1160T98	Balance, 220g x 0.1mg, 110-240V	1/Each	\$2,434.63	\$1,826.00
1160T99	Balance, 310g x 0.1mg, 110-240V	1/Each	\$2,767.39	\$2,075.55

labForce

HDPE Wide Mouth Wash Bottles

labForce High Density Polyethylene (HDPE) Wide Mouth Laboratory Bottles meet the requirements of the most demanding applications. Made from virgin HDPE resin, labForce bottles are durable and break resistant eliminating the shatter risk of glass bottles. The wide mouth opening is perfect for sample collection and allows easy filling and emptying of solids, powders, and viscous liquids. Linerless Polypropylene (PP) caps provide leakproof closure.

- Precision molded for consistent wall thickness and smooth interior and exterior surfaces
- Resistant to most corrosives and suitable for a wide variety of laboratory uses
- Safe for freezer storage down to -100°C
- Plastic laboratory bottles are lighter than glass for improved handling and reduced shipping cost
- Leakproof closure with integrated seal ring



Thomas No.	Description	UOM	List Price	Sale Price
1155L44	Bottle, Round, PP Closure, 30mL	12/Pack	-\$30.90	\$20.10
1155L45	Bottle, Round, PP Closure, 60mL	12/Pack	-\$15.60	\$10.15
1155L46	Bottle, Round, PP Closure, 125mL	12/Pack	-\$18.60	\$12.10
1155L47	Bottle, Round, PP Closure, 250mL	12/Pack	-\$30.95	\$20.15
1155L48	Bottle, Round, PP Closure, 500mL	12/Pack	-\$45.60	\$29.65
1155L51	Bottle, Round, PP Closure, 1000mL	6/Pack	-\$47.00	\$30.55
1155L53	Bottle, Round, PP Closure, 2L	6/Case	-\$148.55	\$96.60
1155L54	Bottle, Round w/handle, PP Closure, 4L	6/Case	-\$210.70	\$137.00
1155L55	Bottle, Square w/handle, PP Closure, 4L	6/Case	-\$240.95	\$156.65

Use promo code
0120
to receive your
special pricing.

labForce

BioRocker™ FX Tube Rocker



- Optimized speed and tilt angle
- 2D "see-saw" motion for thorough, gentle mixing of tubes from 1.5mL to 50mL
- Incubator & cold room safe
- Grooved rubber mat holds tubes in place
- No adapters or mounting hardware

Technical Data

Speed: 24 rpm

Tilt Angle: 30 degrees

Platform Size: 30 x 11cm (12 x 4.5")

Operating Temp. Range: +4°C to 65°C

Weight: 2.2 kg (4.0 lbs)

Manufacturer Warranty: 2 years

Thomas No.	Description	UOM	Price
1154J74	labForce TubeRocker Compact Rocker, 115V	1/Each	\$349.00

labForce

PP Wide Mouth Wash Bottles

labForce Polypropylene (PP) Wide Mouth Laboratory Bottles meet the requirements and expectations of users that value precision. Made from virgin Polypropylene resin, labForce bottles offer excellent chemical resistance and clarity.

- Precision molded for consistent wall thickness and smooth interior and exterior surfaces
- Cap and bottle have semi-buttress threads for better performance than round threads
- USP Class VI resin is FDA -21CFR compliant for food contact
- Autoclavable and approved for use from -40°C to 121°C



Thomas No.	Description	UOM	List Price	Sale Price
1155L56	Bottle, Wide Mouth, Round, PP with PP Closure, 30mL	12/Pack	-\$25.50	\$16.60
1155L58	Bottle, Wide Mouth, Round, PP with PP Closure, 60mL	12/Pack	-\$30.30	\$19.70
1155L59	Bottle, Wide Mouth, Round, PP with PP Closure, 125mL	12/Pack	-\$40.15	\$26.10
1155L60	Bottle, Wide Mouth, Round, PP with PP Closure, 250mL	12/Pack	-\$67.95	\$44.20
1155L61	Bottle, Wide Mouth, Round, PP with PP Closure, 500mL	12/Pack	-\$74.65	\$48.55
1155L62	Bottle, Wide Mouth, Round, PP with PP Closure, 1000mL	6/Pack	-\$68.60	\$44.60

labForce

Cardboard Freezer Boxes with Dividers

labForce cryo/freezer boxes and dividers are made of sturdy white cardboard, with a moisture-repellent coating to protect from the effects of liquid nitrogen or condensation from refrigerated equipment. Boxes are sized to fit standard stainless-steel and polycarbonate storage systems. Small holes on the bottom allow for safe drainage of liquid nitrogen or water. For use in mechanical freezers and liquid nitrogen tanks at temperature to -196°C.

- With a moisture-repellent coating to protect from condensation
- A grid on the inside of the box (81- and 100-place boxes) for easy lining up and location tracking
- Alphanumeric indexing on outer box available



Thomas No.	Description	UOM	List Price	Sale Price
1141Q92	2" Cardboard Freezer Box with 81-Place Cell Divider	10/Pack	-\$20.60	\$12.40
1141Q93	2" Cardboard Freezer Box with 100-Place Cell Divider	10/Pack	-\$20.60	\$12.40
1141Z25	2" Cardboard Freezer Box with 81-Place Cell Divider & Indexing	10/Pack	-\$20.60	\$12.40
1141Z26	2" Cardboard Freezer Box with 100-Place Cell Divider & Indexing	10/Pack	-\$20.60	\$12.40

BUY ONE GET ONE

labForce

Mini Digital Incubators

Receive a FREE labForce Miniature Vortex Mixer when you purchase a Mini Digital Incubator. Offer valid through April 30, 2020.



- Digital temperature control
- Unique air flow design
- Personal sized, 13 x 14.5" footprint
- Accepts bottles & flasks up to 2L
- Two models: heat/cool or heat only
- Optional 3D mini shaker available for sample agitation

FREE labForce Miniature Vortex Mixer!

Thomas No.	Description	UOM	Price
20A00C739	PROMOTION! 1154J59 labForce Mini Incubator (heat only) with FREE 1154J86 Mini Vortexer	1/Each	\$499.00
20A00C740	PROMOTION! 1154J61 labForce Mini Incubator (heat/cool) with FREE 1154J86 Mini Vortexer	1/Each	\$759.00

BUY ONE GET ONE

labForce

Gradient Thermal Cycler

Receive a FREE labForce Mini Dry Bath when you purchase a Thermal Cycler. Offer valid through April 30, 2020.



- Intuitive touch screen control
- Program Wizard- quick, easy & foolproof
- Unique multi-format block
 - 0.2mL tubes or strips (96)
 - 0.5ml tubes (39)
 - 96 well plate (1)

FREE labForce Mini Dry Bath!

Thomas No.	Description	UOM	Price
20A00C741	PROMOTION! 1154J63 labForce Microcentrifuge with FREE 1154J48 Mini Dry Bath	1/Each	\$4,999.00



Disclaimer: Pricing and offers subject to change without notice.

Thomas
Scientific

ThomasSci.com
800.345.2100

Thomas Scientific, LLC Family of Companies

DENVILLE
SCIENTIFIC

PHENIX
RESEARCH PRODUCTS

ESK SCIENTIFIC

AMERICAN
CLEANSTAT