

FREE Equipment

with purchase of item below*
(Free with items: 1227U62, 1165T06, 1165T08, 1165T10, 21A00H426, 1227W13/1227W14, 20A00R247, or 1159P17)



Value: \$414.49

Clik-n-Lok™
rotor attachment



Value: \$410.63

MicroSterilizer



Value:
\$429.43

Tube Rotator
Accepts tubes 1.5 to 50ml

Your choice of:

MyFuge™ Mini
Item: 1222U52

BactiZapper™
Item: 22A00B664

Roto-Mini™
Item: 1159P32

BENCHROCKER™ 2D Variable Speed 2-D Rocker

\$981

- Adjustable speed and tilt
- Ideal for staining and general mixing
- Safe for incubator or cold room use

Speed: 2 to 30 rpm
Platform: 14 x 12 in.
Tilt Angle: 0 to ±30°



Optional Stacking Platform

1227U62	BenchRocker™ 2D Variable Speed Rocker	\$981.44
1227U67	Optional stacking platform	\$162.44
1227U69	LaBungee™ elastic tie downs with hooks, 4/pk.	\$37.89

SURETEMP™ Incubators with SureCheck™ Data Logging

- Precision temp. control up to 75°C
- Built in, SureCheck™ data logger
- Dual convection, gravity or forced air
- Stainless steel interior

Chamber Dimensions:
40L: 13.7 x 13.7 x 13.7 in.
70L: 17.7 x 13.7 x 17.7 in.
130L: 21.7 x 17.7 x 21.7 in.



1165T06	SureTemp Incubator, 40L	\$1,871.77
1165T08	SureTemp Incubator, 70L	\$2,412.90
1165T10	SureTemp Incubator, 130L	\$3,217.59

MULTI-THERM™ COOL - HEAT - SHAKE

- Heats, cools and shakes
- Variety of blocks for 0.2 to 50ml tubes
- Digital programming with sequence link function *

*Only available in 1227W13 model

Temp Range: Ambient -20° to 100°C*
Speed Range: 200 to 1500rpm
*1227W14 model, minimum temp is ambient +5°C



1227W13	MultiTherm™ Shaker, heat and cool	\$3,446.68
1227W14	MultiTherm™ Shaker, heat only	\$2,619.87
1227W18	Block, 35 x 1.5ml	\$368.33
1227W19	Block, 35 x 2.0ml	\$368.33
1227W15	Block, 96 x 0.2ml or 1 x PCR plate	\$438.23
1227W21	Block, 12 x 15ml	\$519.53
1227W22	Block, 6 x 50ml	\$519.53

ROTO-THERM™ INCUBATED TUBE ROTATOR

- Temperature control up to 60°C
- Includes tube holders for 0.5 to 50ml
- Easy, in-lab temperature calibration
- Available with variable speed & mixing modes
 - Oscillating
 - Pausing
 - Rocking

Included Holders

- 11 x 1.5-2ml (2/pk.)
- 5 x 5-7ml (2/pk.)
- 4 x 10-15ml (2/pk.)
- 6 x 50ml* (1 ea.)

*Also accepts 24 x 0.5ml

Temp Range: Ambient +5° to 60°C
Speed Range: 5 to 70 rpm
Max. Capacity: 6 x 50ml, 16 x 15ml, 24 x 0.5ml, 44 x 1.5/2.0ml, 96 x 0.2ml



1159P17	RotoTherm Plus w/included tube holders (variable speed)	\$803.49
1159P33	Optional rocking platform tube holder for 96 x 0.2ml (8 x PCR strips)	\$88.41

BenchPress

Vol. 1.1 2022
Promo: BP1.22

**Thomas
Scientific**

Shakers Incubators Centrifuges Balances Rotators Vortexers DryBaths Rockers and MORE!

BeadBlaster™ 24 Homogenizer Series

- Completely lyse, grind & homogenize a variety of sample types
- Available with REFRIGERATION to protect delicate samples
 - True refrigeration - no dry ice or LN₂ needed
 - Space saving design
- Process up to 24 samples simultaneously
- Wide variety of grinding media in prefilled 2ml & 5ml tubes

Specifications (1201P15 and 20A00J667)

Capacity: 24 x 2ml tubes
Speed range: 4.00 - 7.00 m/sec
D2400-R Only:
Cooling: Compressor
Temp. Range: -10°C to Ambient

Special Offer

1201P15	BeadBlaster 24 Homogenizer	\$8,875.65	\$11,999.00
20A00J667	BeadBlaster 24R with refrigeration	\$13,954.58	
22A00B643	Optional tube holder for 5ml tubes	\$413.11	
1211U65	Prefilled 2ml tubes, 0.5mm zirconium, 50/pk.	\$196.75	
1211U67	Prefilled 2ml tubes, 1.5mm zirconium, 50/pk.	\$197.74	
1211U68	Prefilled 2ml tubes, 3mm zirconium, 50/pk.	\$214.22	
1211U70	Prefilled 2ml tubes, 2.8mm steel, 50/pk.	\$197.74	
21A00K863	Prefilled 5ml tubes, 3mm zirconium, 50/pk.	\$219.00	
1159P01	Prefilled 5ml tubes, 2.8mm steel, 50/pk.	\$233.88	

Additional grinding media is available in prefilled tubes and loose bulk. Ask for details.

Also Available:

1211U60
3 Position

1165T60
6 Position



QuadCount™ Automated Cell Counter

- Reduce counting errors, save time and improve efficiency
- Intuitive software and touch-screen interface
- Convenient and cost-effective QuadSlides™
- Quick, Regular and Precise count modes (up to 3.6µl count volume)

21A00B769	QuadCount™	\$3,879.75
21A00B771	QuadSlide™, 50/pk.	\$156.45

Replaces manual cell counting...ask for a demo today!

*Promo Code BP1.22

BioClave™ Autoclave Series

- Designed specifically for the sterilization of laboratory research tools
- Fully automatic operation - just press start
- Built in water tank - no external connection needed
- 8L, 18L and 28L models

Accepts bottles and flasks up to 2L



1227W09	8L BioClave	\$3,522.69	21A00Q414	28L BioClave	\$11,545.20
21A00Q412	18L BioClave	\$6,240.25	1142K09	28L BioClave (220V)	\$11,776.11

SureTherm™ CO₂ Incubators with INCUVIEW™ LCI

NEW


20A00V149

20A00V148



*Promo Code BP1.22

- Internal microscope option for **Live Cell Imaging**
- Precision temperature control with 6-sided heating
- Available with:
 - Digital CO₂ and O₂ control
 - High heat decontamination
 - Split inner door



Capture time-lapse images

Special Offer

20A00V140	ST-45 CO ₂ Incubator with two shelves	\$5,664.75	\$5,295
20A00V142	ST-180 CO ₂ Incubator with three shelves	\$6,415.41	\$6,295
20A00V144	ST-180 CO ₂ Incubator w/high heat decontamination	\$8,127.41	\$7,895
20A00V146	ST-180 CO ₂ Incubator w/high heat decontamination split window door & O ₂ control	\$11,123.41	\$10,595
20A00V148	IncuView LCI Integrated microscope shelf (14 x 11.5 in.)	\$1,690.20	
20A00V149	IncuView LCI Integrated microscope shelf (20 x 17 in.)	\$1,776.60	

NEW! View cells without disturbing their CO₂ environment...
On the front display, or remotely from home!

**Thomas
Scientific**



ThomasSci.com
833.544.SHIP (7447)
CustomerService@thomassci.com

Connect With Us:



*Promo Code BP1.22 Valid Through 06.30.22. Pricing, availability, and specifications subject to change without notice.



BEADBATH Duo™

THE ONLY BATHS DESIGNED FOR BEADS OR WATER

- Two separate modes for BEADS or WATER
- Quik-Cal™, simple in-lab calibration function
- Offers all of the advantages of a "Bead Bath"
 - Eliminates constant refills (burn-out protection)
 - Easy (accessory free) securing of tubes and bottles
 - Allows angled incubations
 - Runs constant (no wait for heat-up)

20A00U450 BeadBath Duo™ 2L with lift-off lid \$744.12
 20A00U452 BeadBath Duo™ 10L with hinged lid \$1,269.00
FREE UV Clave™ with purchase of any BeadBath Duo



Choose Bead or Water Mode *Promo Code BP1.22

Special Offer! Item: 20A00R635
FREE: UV CLAVE™ Ultraviolet Chamber



Value:
\$511.21

Inactivates DNA, RNA, viruses & bacteria on PCR vessels

SMARTREADER™ 96

Your go-to instrument for ELISA and absorbance assays



- Intuitive, 7 inch, color touch screen control panel
- Filter based, wavelength range from 340-750nm (405, 450, 492, 630nm included)
- Operates as a stand alone system, with USB flash drive for data transfer
- Integrated, multi-speed plate shaker

Also Available:

20A00R629
 Microplate Washer



1165Q91 SmartReader 96 Microplate Absorbance Reader \$6,356.39
 20A00L661 SmartReader 96 Microplate Absorbance Reader with Incubation \$6,440.50
 1165Q93 SmartReader Thermal Printer \$358.33

Roto-Mini™ Plus

- Includes tube holders for common tube sizes
- Digital speed control
- Advanced shaking/mixing modes:
 - Mix (oscillating)
 - Pause
 - Rock

\$375.39



1159P34 Roto-Mini™ Plus Variable Speed Rotator \$375.39



MC-24™ Touch

High Speed Microcentrifuge

- Full color touch screen
- mySpin™ 1 button program recall
- COMBI-Rotor™ accepts tubes and strips
- Compact, only 9 x 12 x 7.75 in.
- Optional rotors for 5ml tubes and 2ml cryovials



Touch Control

Max. Speed: 13,500rpm / 16,800xg
 Max Capacity: 24 x 1.5/2.0ml
 2 x PCR Strips (16 x 0.2ml)



\$2,375.11

Special Offer! Item: 1220P38
FREE: myFuge™ 12 Mini Centrifuge

Value:
\$331.23



*Promo Code BP1.22

1160X90 MC-24 Touch Microcentrifuge with 24 place COMBI-Rotor \$2,375.11
 21A00K866 Optional 8 x 5ml rotor for MC-24 Touch \$1,004.98
 21A00K867 Cryovial adapter for C2417-R5 rotor, 8/pk \$125.35

FREE myFuge™ 12 with purchase of MC-24 Touch Microcentrifuge

HERMLE Benchmark



As low as \$5,800

Z216-MK High Speed Microcentrifuge

REFRIGERATED

- Temp Range -20°C to +40°C
- High speed, up to **21,631xg**
- High capacity, up to **44 x 1.5/2.0ml**
- **10 rotor options for tubes up to 5ml**
- Innovative **COMBI-rotor** holds:
 24 x 1.5/2.0ml tubes **and** 4 x PCR strips!
- **5 year warranty!**



44 x 1.5/2.0ml



24 x 1.5/2.0ml
 and 4-PCR strips

BUNDLE ORDERING INFORMATION

Reg. **Bundle**

1159P62	Z216-MK Refrigerated w/ 44 place rotor	\$7,012	\$5,800
1144F72	Z216-MK Refrigerated w/ COMBI-Rotor™	\$7,714	\$6,200

*Ambient versions also available

INCUSHAKER™

SHAKING INCUBATORS

- **MAGic™** Clamp accessories for instant-magnetic clamp exchange
- Available with CO₂ control or cooling
- **NEW!** Adjustable clamp for 125, 250 and 500ml flasks (4 FREE with purchase)

Limited Time Offer*
FREE!

4ea. **MAGic Clamp (Adjustable)**


PATENTED Item: 1177X57
 (Accepts 125ml, 250ml and 500ml Flasks)

*Promo Code BP1.22

INCUSHAKER Mini

Speed Range: 30 to 300rpm
 Temp Range: Amb. +5° to 70°C
 Platform Dims: 9.5 x 11.5 in.
 Universal Platform: 1222U14
 Flask Capacity: 2 x 1L (1222U28)
 5 x 500ml (1222U15)
 8 x 250ml (1222U16)
 12 x 125ml (1222U17)

\$2,587.35



Item: 1156Z77

INCUSHAKER 10L & 10LR

Speed Range: 30 to 300rpm
 Temp Range: (1203F37, Cooling Model)
 Amb. +5° to 60°C
 Amb. -15° to 60°C
 Platform Dims: 17.8 x 17.8 in.
 Universal Platform: 1207X19
 Flask Capacity: 5 x 2L (1222U27)
 9 x 1L (1222U28)
 16 x 500ml (1222U15)
 25 x 250ml (1222U16)
 36 x 125ml (1222U17)

\$8,029.50 / \$11,295.38



Item: 1213F79 & 1203F37 (Cooling)

INCUSHAKER Mini CO₂

Speed Range: 30 to 300rpm
 Temp Range: Amb. +5° to 60°C
 CO₂ Range: 0 to 20%
 Platform Dims: 9.5 x 11.5 in.
 Universal Platform: 1222U14
 Flask Capacity: 5 x 500ml (1222U15)
 8 x 250ml (1222U16)
 12 x 125ml (1222U17)

\$10,449.71



Item: 1177X58

my BLOCK™ Digital Dry Baths

- Digital temperature control up to 105°C
- Hinged lid improves uniformity & conserves energy
- External temperature probe (optional)
- Available with **Quick-Flip™** universal blocks for microtubes (0.2 to 2.0ml), PCR strips & PCR plates

1203F25	myBlock I single chamber dry bath	\$580.31
1203F27	myBlock II dual chamber dry bath	\$722.47
1203F29	myBlock I single chamber dry bath, with 1x BSWCMB	\$745.48
1203F31	myBlock II dual chamber dry bath, with 2x BSWCMB	\$1,046.83
1203F35	Optional external temp. probe, for sample temp. control	\$246.87

(Blocks sold separately)

Block Ordering Information

1203F36	Quick-Flip, 24 x 1.5/2.0ml, 14 x 0.5ml and 6 x PCR Strips / 48 x 0.2ml	\$210.19
1213L24	48 x 0.2ml (6 x PCR Strips)	\$170.89
1227U84	24 x 1.5ml	\$170.89
1227U85	24 x 1.5/2.0ml	\$170.89
1227U89	20 x 13mm	\$170.89
1227U87	12 x 15ml conical	\$170.89
1227U90	5 x 50ml conical	\$170.89

(Additional blocks available, please visit our website for a complete listing)

QuickFlip™ Reversible Block

1203F36



24 x 1.5/2.0ml & 14 x 0.5ml
 48 x 0.2ml or 6 x PCR Strips*

*Place 2 blocks together in a dual dry bath to accept a 96 well PCR plate

