

Glass Disposal Boxes

Floor and Bench Models *Safely handle broken glass and spilled liquids!*

- Sturdy corrugated cardboard
- 2 mil thick polypropylene liner
- Two models: bench and floor

Keep the lab safe with Thomas glass disposal boxes. Constructed from rugged corrugated cardboard and supplied with a 2 mil thick puncture-resistant polyethylene liner, the unit safely contains broken glass and spilled liquids. When full, simply pull the safety cap into place and dispose of the entire unit. Striking graphics make the glass disposal box easy to find. Box comes in two versions: Benchtop (8" x 8" x 10") and Floor (12" x 12" x 27") which hold 5.5 and 40 lbs respectively.



Thomas No.	Description	Pack
1236E00	Glass Disposal Box, Floor	6
1236E01	Glass Disposal Box, Bench	6

Glass Disposal Box Cover

- Dishwasher-safe and autoclave-safe at 112°C (250°F)
- Reusable cover opens quickly and easily by pressing the built-in release tab
- Fits any disposal carton that has a 12" x 12" opening

Provides a wider opening (8.25 x 9.5") than cardboard lids, permitting safer and more sanitary deposit of wastes. at 112°C (250°F).



Thomas No.	Description	Each	Case
1225Z68	Glass Disposal Box Cover	1	8