

Thomas
BRAND

Portable Density System

This self-contained kit includes milligram balance and a complete solids and liquids density set, with thermal printer.

Balance model S303, capacity 310g, resolution 0.001g

- Cylindrical glass draftshield Ø135mm high 70mm
- Weighing pan Ø 80mm
- Internal rechargeable battery
- Switching power supply 9V 1A with interchangeable plugs USA-EU-UK
- Weighing pan Ø 80mm
- Calibration weight 200g F1
- Software for automatic calculation of solid and liquid density



Density set for solids and liquids

- Measuring glass cup 150mm high Ø 50mm 250ml
- Measuring glass cup 95mm high Ø 65 mm 250ml
- Thermometer
- Glass sinker
- Sample holder for floating/non floating samples
- Hanging arm for sample holder or sinker
- Sample holder/sinker support
- Measuring glass cup support

Thermal printer DPP

- Internal rechargeable battery
- Switching power supply
- RS232 cable
- Thermal paper roll

| Thomas No. | Description | Qty. /Pk |
|----------------|--------------------------------------------|----------|
| 1200J68 | Density Determination Kit, Portable System | 1 |
| 1200J69 | DPP-250 Printer with cable | 1 |
| 1200J70 | Roll of Paper for DPP-250 Printer | 1 |
| 1200J71 | Test Dishes for Moisture | 80 |