

25% OFF Heavyweight Cleanroom Wipes



Designed for critical cleaning applications, M1 wipes are made from 100% continuous filament polyester in a durable double-knit construction. Manufactured, washed, and packaged in-house for exceptional cleanliness and quality control.

Composition & Attributes

- 100% Hi Strength Polyester fiber (continuous filament, double knit)
- Ultra Low particle and fiber generation
- Laser Sealed Edges for fiber & particle retention
- Free of Silicon, Amides, and DOP contamination
- Solvent safe double bag cleanroom packaging
- Resists abrasion when used with rough surfaces
- Works well with IPA and other cleaning solvents

Use promo code **0726 at checkout to receive special promotional pricing.**

Promotion expires 8/31/2026. Disclaimer: Pricing and offers subject to change without notice.

Thomas No.	Description	Qty.	Price	New Price
21A00P823	M1-3400LSSE Cleanroom Wipes, Heavyweight Polyester Woven Laser Sealed Edge, Size 12 x 12	100/PK	72.19	54.14
		1000/CS	721.90	541.40



customerservice@thomassci.com



thomassci.com

