



Pre-Saturated Alcohol Wipes

PC-SAT

- Contains 160 - 6" x 7" Wipes
- Pre-Saturated with 70% Isopropyl Alcohol / 30% DI Water
- Polyester / Cellulose Non-Woven Material
- Foil Sealed for Saturation Assurance
- Re-Sealable Closure for Freshness

Thomas No.
20A00M182

Thomas Scientific, LLC
1654 High Hill Road • Swedesboro, NJ 08085
ThomasSci.com

General Customer Service: 833.544.7447 | Production Customer Service: 833.750.1234
Made in China

Isopropyl Alcohol/Alcohol Isopropílico



Warning:

Wipers contain flammable liquid and vapor, may cause respiratory, skin irritation, and eye injury. Harmful if inhaled or swallowed.

Precautions:

Keep away from heat, sparks, and open flame. Use only in well-ventilated, smoke free areas. Keep away from sources of ignition. Avoid contact with skin and eyes. Wear proper protective clothing, gloves, and eye protection.

If contact on skin or contact with eyes: rinse immediately with plenty of water and seek medical advice.

If inhaled, move to fresh air. If not breathing, provide artificial respiration and get medical attention.

Emergency Phone: **CHEMTREC, U.S.: (800) 424-9300**
International: (703) 527-3687

Read the material safety data sheet before using this product.
For industrial use only.

Advertencia:

Paños con líquido inflamable y vapor. Puede causar dificultad respiratoria e irritación en ojos y piel. Nocivo si se inhala o ingiere.

Precauciones:

Manténgase alejado de fuentes de calor, chispas o llamas directas. Utilice solo en zonas ventiladas y áreas libres de fumadores. Mantenga el producto alejado de fuentes de ignición. Evite el contacto con la piel y ojos. Utilice indumentaria, guantes y protección para ojos adecuada.

En caso de contacto con ojos y piel: lavar con abundante agua el área afectada y obtenga atención médica.

En caso de inhalación: traslade a la persona afectada al aire libre. Si no respira o presenta dificultad, suministrar respiración artificial y obtenga asistencia médica.

Teléfono de emergencia:
CHEMTREC, U.S.: (800) 424-9300
Internacional: (703) 527-3687

Leer hoja de seguridad antes de utilizar este producto.
ÚNICAMENTE PARA USO INDUSTRIAL.