

Customer Change Notification

Thomas Scientific is Migrating ThermoGrid™ PCR Tubes to NEW Amplifyt™

Thermo Grid™



Amplifyt™
powered by labForce

Thomas Scientific is proud to launch an expansion of products offered through our own "labForce" brand. As always, you can trust that we took great care in creating the best products to support our customers in their research, while also maintaining the best possible pricing. Thomas Scientific stands behind our new Amplifyt™ products and ensure they will not only meet your research requirements but exceed your expectations.

Amplifyt™ PCR tubes are *equivalent* in fit and function to the original ThermoGrid™ tubes while also providing these *additional superior benefits*:

EASY Handling:

- Amplifyt™ tubes feature rigid linkages



IDEAL for qPCR :

- Amplifyt™ tubes offer flat, optically-clear caps



INCREASED Clarity:

- Ultra-clear for easy visualization

EASY to Break:

- Single point linkage
- Snap off only the tubes you need

IMPORTANT: While ThermoGrid™ and Amplifyt™ are equivalent in function, caps and tubes from each brand are not compatible or interchangeable. Please use Amplifyt™ caps with Amplifyt™ tubes to ensure optimal sealing.

12/06/2019

Customer Change Notification

Thomas Scientific is Migrating ThermoGrid™ PCR Tubes to NEW Amplifyt™

ThermoGrid™



Amplifyt™
powered by labForce

What you can expect during the transition from ThermoGrid™ to Amplifyt™:

- This transition will occur over the next several months. Each part number will transition at a different time – depending on stock.
- Below you will see the old ThermoGrid™ parts and the new Thomas Scientific part number for the equivalent Amplifyt™ product will receive in the future.
- The switch will occur automatically when our stock of any part number is depleted. You will be notified of required changes to your purchase order.



Individual PCR Tubes

ThermoGrid™ Thomas No.	ThermoGrid™ Qty per Pack	NEW Amplifyt™ Thomas No.	Amplifyt™ Qty per Pack
1158R16	1,000	1148A00	1,000
1158U19	1,000	1149Z98	1,000
1230Z96	1,000	1148A00	1,000

PCR Strip Tubes with Attached Caps

ThermoGrid™ Thomas No.	ThermoGrid™ Qty per Pack	NEW Amplifyt™ Thomas No.	Amplifyt™ Qty per Pack
1159N19 ¹	125	1148A28 ¹	120
1148A23	120	1148A24	120
1148A27	120	1148A28	120
1156C74	120	1148A20	120
1159N20 ²	125	1148A16 ²	120
1159V37	120	1148A16	120
1159V38	120	1148A25	120

¹Amplifyt™ caps on 1148A28 are optically clear and flat for qPCR applications

²Amplifyt™ strip is easy to break (compared to ThermoGrid™ rigid strip)

PCR 12-Strip Tubes with Strip Cap Packs

ThermoGrid™ Thomas No.	ThermoGrid™ Qty per Pack	NEW Amplifyt™ Thomas No.	Amplifyt™ Qty per Pack
1149Z77	125	1149K07	125
1158U18	125	1149K07	125

PCR 4-Strip Roto-Gene™-Style Tubes

ThermoGrid™ Thomas No.	ThermoGrid™ Qty per Pack	NEW Amplifyt™ Thomas No.	Amplifyt™ Qty per Pack
1148G56	250	1139Y44	250

PCR 8-Strip Tubes with Strip Cap Packs

ThermoGrid™ Thomas No.	ThermoGrid™ Qty per Pack	NEW Amplifyt™ Thomas No.	Amplifyt™ Qty per Pack
1148G52	125	1149Z65	125
1149Z55	125	1149Z65	125
1125Z57	125	1149Z65	125
1149Z58	125	1139Y39	125
1158U17	120	1149Z65	125

Separate PCR Strip Caps and Strip Tubes

ThermoGrid™ Thomas No.	ThermoGrid™ Qty per Pack	NEW Amplifyt™ Thomas No.	Amplifyt™ Qty per Pack
1158R13	125	1149Z13	125
1158T85	120	1149Z13	125
1158U47	120	1149Z34	125
1148G51	125	1149Z36	125
1158R15	125	1149Z36	125
1219F17	125	1149Z36	125

*Highlighted item indicates change in quantity when moving to Amplifyt™

¹While ThermoGrid™ and Amplifyt™ are equivalent in function, caps and tubes from each brand are not compatible or interchangeable. Please use Amplifyt™ caps with Amplifyt™ tubes to ensure optimal sealing.

