

Petri Dishes

***Thomas Disposable Sterile Petri Dishes
are available in stackable and slippable forms***

- Virgin biomedical grade polystyrene
- Optically clear
- Vented lids for air circulation
- Radiation sterilized



Stackable Thomas Petri Dishes are ideal for general laboratory use while Slippable Thomas Petri Dishes are designed for automated filling systems. Vented lids permit air circulation during incubation to reduce condensation build-up. Petri dishes are optically clear for easy viewing.

Dishes are available in the following nominal sizes: 60 x 15mm, 100 x 10mm, 100 x 15mm, and 150 x 15mm.

Thomas No.	Description	Case
1188N80	Petri Dish, 100 x 15mm, Slippable, Sterile, 500	500
1188N81	Petri Dish, 100 x 15mm, Stackable, Sterile, Bulk, 25/500	100
1188N82	Petri Dish, 60 x 15mm, Stackable, Sterile, 25/500	500
1188N83	Petri Dish, 150 x 15mm, Slippable, Sterile, Bulk, 10/100	500