

# Serological Pipets

## Individually Wrapped, Paper/Plastic



- Available with 7 capacities of 1.0, 2.0, 5.0, 10.0, 25.0 and 50.0ml
- 1.0, 2.0, 5.0, 10.0ml are stretched, while 10.0, 25.0 and 50.0ml are ultrasonically welded at tip and mouth-piece
- Graduations are calibrated for accurate dispensing to within  $\pm 1\%$
- Different color ring for easy identification
- Bidirectional graduations on the pipets provide added applicability
- Negative graduation allows additional working volume
- All pipets are supplied with a filter plug
- Strict consistency tested
- Non-pyrogenic



Thomas No.	Mfr. No.	Description	Pack/Case Qty.
<b>21A00B813</b>	GSP010001	Serological Pipet, 1mL, Individually Wrapped, Yellow Band, Sterile	500 Case
<b>21A00B814</b>	GSP010002	Serological Pipet, 2mL, Individually Wrapped, Green Band, Sterile	500 Case
<b>21A00B815</b>	GSP010005	Serological Pipet, 5mL, Individually Wrapped, Blue Band, Sterile	200 Case
<b>21A00B816</b>	GSP010010	Serological Pipet, 10mL, Individually Wrapped, Orange Band, Sterile	200 Case
<b>21A00B817</b>	GSP010025	Serological Pipet, 25mL, Individually Wrapped, Red Band, Sterile	150 Case
<b>21A00B818</b>	GSP010050	Serological Pipet, 50mL, Individually Wrapped, Purple Band, Sterile	100 Case