



# Labels for the Cannabis Industry



GA International provides cannabis labels and labeling solutions for all your cannabis industry needs. We offer plant labeling systems, tube, vial, and container labels for QC testing, as well as custom branding and packaging labels.



## Packaging Labels for Cannabis Brands, Distributors, & Retailers

- Full-color custom printing of logos and designs made to your specifications, in digital, flexographic, and thermal formats
- Shorten turn around times, reduce wait times, and be sure that your product will be ready to ship ASAP
- Wide selection of labels that conform to a variety of surfaces, including plastic bags, jars, and bottles, all of which may be used to package dried flowers, seeds, edibles, tinctures, oils, or pre-rolled products
- Our labels meet FDA regulations for indirect food contact and are RoHS compliant
- State and country-specific warning and compliance labels are available
- Tamper-evident labels can be used on bags and jars to provide additional security and to ensure that they remain in compliance with both local and federal regulations





## Plant Identification Labels for Cannabis Growers

- Identifies the plant, not the pot, ensuring every plant remains clearly labeled from seed to sale, even if repotted or moved
- Durable and waterproof, the labels will remain attached and legible, even in damp, humid environments
- Can be printed with various high-quality barcodes and variable information for enhanced tracking and traceability
- Available in a variety of colors, perfect for color coding and breeding of unique plant strains

SKU#	Inches (W x H)	mm (W x H)	Core	Labels/Roll
B-OL-1C1-1	2.5" x 1.25"	63.5 x 31.8	1	1,000
B-TTOL-33C1-1	3.063" x 1.25"	77.8 x 31.8	1	1,000



## Tube, Vial, & Dish Labels for Testing Labs

- We provide durable labels that will adhere to small diameter tubes and vials of different sizes, including small LC vials, 50 ml Falcon tubes, and culture plates
- These labels resist spraying and wiping with alcohols, such as ethanol and isopropanol, as well as other chemicals and detergents
- The labels will also withstand deep-freeze conditions if cold storage is required (e.g, for batch retention)
- Thermal-transfer printable, allows the labels to be printed with barcodes



SKU#	Inches (W x H)	mm (W x H)	Core	Labels/Roll	Surface
FJT-7C1-1	1" x 0.5"	25.4 x 12.7	1	1,000	Vial
FJT-29C1-1	1" x 1"	25.4 x 25.4	1	1,000	Vial
FJT-57C1-1	1.125" x 0.75"	28.6 x 19.1	1	1,000	Vial
FJT-97C1-1	1.28" x 0.625"	32.5 x 15.9	1	1,000	Vial
FJT-28C1-2	2" x 1"	50.8 x 25.4	1	2,000	Tube
FJT-169C1-4	1.25" x 0.25"	31.8 x 6.4	1	4,000	Plate
FJT-11C1-2	2.64" x 0.276"	67.1 x 7	1	2,000	Plate



## Printers, Scanners, & Software

Equip your grow-op, lab, or business with a complete labeling solution from GA International. GA offers printers, ribbons, barcoding software, scanners, and everything else you'll need to securely and accurately identify your plants, samples, and packages. Ready-made printing kits are available, or you can build your own cannabis printing solution with the help of our knowledgeable scientific experts.

- Contact your local Thomas Scientific Sales Representative to assist you with more sizes, colors, core sizes, configurations, and matching ink ribbons.