

Safety SPOTLIGHT

**Thomas
Scientific**

We Believe You Are Important
— How Can We Help?

thomassci.com

global@thomassci.com

Inside US **800.524.0018**
Outside US **856.467.2000**

Think **Thomas** for All
of Your **Safety Needs!**



Summer Edition 2017

Use **Promo Code 0217** to receive the featured pricing. Valid until 12/31/17.

**Thomas
Scientific**

Workrite[®]
FR

Workrite[®]
FR CP

THE RIGHT FIT FOR EVERY LAB

Innovations in Laboratory Safety



Workrite FR offers different styles of lab coats to meet various needs. All are flame-resistant (FR), but there are other key differences. Newer styles are made with a lighter-weight fabric, making them more breathable and comfortable. They also have knit cuffs, which were added based on feedback from the University of California. These cuffs keep foreign material out and prevent accidental spills.

The latest innovation is the Workrite FR/CP Lab Coat, which combines the FR properties of DuPont Nomex[®] IIIA with the chemical-splash protection (CP) of Westex ShieldCXP[™]. The first of its kind, this lab coat simultaneously protects against short-term thermal hazards and inadvertent liquid chemical splashes while also providing superior comfort and flexibility.

Chemicals that ShieldCXP is generally effective against include aqueous chemicals and mixtures, and polar organic solvents. For more information, refer to workritefrcp.com

Featuring fabrics from:

DuPont[™]
Nomex[®]

Westex
ShieldCXP[™]

© Workrite Uniform Company, Inc.
Workrite Uniform is an affiliate of Williamson-Dickie Manufacturing Company
Nomex[®] is a registered trademark of DuPont. | ShieldCXP[™] is a trademark of Westex, Inc.

LAB COAT
Thomas Nos.
1169L19 – 1169L28



**LAB COAT
W/KNIT CUFFS**
Thomas Nos.
1169L39 – 1169L48



**FR/CP
LAB COAT**
Thomas Nos.
1169L29 – 1169L38



| | | | |
|--|---|---|---|
| Flame-resistant (FR) protection | ✓ | ✓ | ✓ |
| Styles for men and women | ✓ | ✓ | ✓ |
| Concealed snap front closure | ✓ | ✓ | ✓ |
| Pass-through access to inner clothing | | ✓ | ✓ |
| Knit cuffs for improved comfort and protection | | ✓ | ✓ |
| Lightweight 4.5-ounce DuPont Nomex [®] IIIA fabric | | ✓ | ✓ |
| Chemical-splash protection (CP) with Westex ShieldCXP [™] | | | ✓ |

Inside US **800.524.0018** Outside US **856.467.2000** | thomassci.com | global@thomassci.com



Think Thomas for All of Your Safety Needs!

Thomas
Scientific

PIG New Pig

PIG® Spill Kit in Bucket

- Polyethylene bucket with easy-open, screw-on lid can be used to ship or temporarily store used absorbents
- PIG BLUE Socks stop spreading spills; PIG Mat absorbs quickly
- Prepacked container provides easy access to contents
- Large, color-coded label makes identification fast and accurate
- Tamperproof seals help prevent pilfering of spill response supplies
- UN Rating 1H2/Y30/S for Packing Groups II and III
- Use as refill for wall-mount holder (sold separately)
- Only PIG Spill Kits feature PIG Absorbents proven for rapid response.
- For information on custom spill kits, just call 1-800-HOT-HOGS (468-4647)



| Thomas No. | Description | Each | List |
|------------|-----------------------------------|------|-------|
| 1205U97 | PIG® Spill Kit in Bucket | 1 | 73.33 |
| 1205U98 | PIG® Oil-Only Spill Kit in Bucket | 1 | 85.99 |

Thomas
Scientific


Universal Spill Kit

- Neutralizes acids, bases organics, hydrocarbons, water reactive liquids, poisons, air reactive liquids, oxidizers, reducing agents and a host of other types of chemicals
- Color-coded and designed to change colors
- Non-toxic
- Safe to handle
- Easy-to-store
- Easy-to-use
- Rust inhibitor prevents degradation of hazmat shipping and storage containers



Ampho-Mag™ is a rapid response powder that safely adsorbs and treats any spill and produces less heat and spattering than any other neutralizer available. A gas neutralizer and adsorbent, it prevents pressure build-up in containers.

| Thomas No. | Description | Each | List | Sale |
|------------|---------------------|------|--------|--------|
| 8238N05 | Universal Spill Kit | 1 | 193.67 | 141.92 |



Looking for signage?

Thomas Scientific offers a wide variety of signs from National Marker (NMC) to keep your workplace safe.

Custom signage also available!

Ask your Thomas Scientific Sales Representative!

Visit thomassci.com
Search National Marker



JUSTRITE

Sure-Grip® EX Under-counter Corrosives/Acid Steel Safety Cabinets

- Safely store corrosives under a workbench or lab counter for quick and easy access
- These cabinets have the same quality safety features as the cabinets for flammables
- Double-wall construction, dual vents, grounding wire connections, adjustable shelves, leakproof sills, three-point self-latching doors and leveling feet
- Manual-close door(s) opens to a full 180 degrees and self-latches when pushed closed
- Superior puncture & abrasion resistance



| Thomas No. | Description | Each | List | Sale |
|------------|---|------|---------|--------|
| 8246B12 | Sure-Grip EX Undercounter Corrosives/Acid Steel Safety Cabinet, Blue, 1 shelf, 2 manual-close doors | 1 | 1081.39 | 869.00 |

Sure-Grip® EX Corrosives/Acid Steel Safety Cabinets

Research environments using many different types of chemicals have special storage requirements. Whether using acids, bases or solvents, it is critical that incompatible liquids are segregated to avoid disastrous explosions or fire. The SDS (Safety Data Sheet) should be reviewed to determine chemical characteristics and recommended storage practices.

- Blue steel safety cabinets for corrosives have the same quality safety features as the cabinets for flammables
- Manual-close door(s) opens to a full 180 degrees and self-latches when pushed closed
- Include polyethylene trays attached to galvanized steel shelves and a separate polyethylene liner for the bottom sump resist aggressive chemicals
- Self-close door(s) shuts and latches automatically when a fusible link melts at 165°F (74°C) under fire conditions
- Manual close models meet OSHA and NFPA 30




| Thomas No. | Description | Each | List | Sale |
|------------|---|------|---------|--------|
| 8246B16 | Sure-Grip EX Undercounter Corrosives/Acid Steel Safety Cabinet, Blue, 22 gal, 1 shelf, 2 manual-close doors | 1 | 1131.49 | 909.27 |



Polyethylene Safety Disposal Cans

Eagle Poly Disposal Safety Cans provide an intermediate station for temporary storage of liquid chemical wastes.

- Cans are FM approved and feature a wide 4" spout
- Safety cans have a cap that locks open to 80 degree angle for easy filling
- High density polyethylene (HDPE) construction resists dents and punctures
- Stainless steel fittings and flame arrestor, ideal for aggressive chemicals and corrosives



| Thomas No. | Description | Each | List | Sale |
|------------|--|------|--------|--------|
| 1197B16 | Disposal Can, Polyethylene, 3 Gallon, Red | 1 | 261.14 | 167.88 |
| 1197B17 | Disposal Can, Polyethylene, 3 Gallon, White | 1 | 261.14 | 167.88 |
| 1197B18 | Disposal Can, Polyethylene, 3 Gallon, Yellow | 1 | 261.14 | 167.88 |
| 1197B19 | Disposal Can, Polyethylene, 5 Gallon, Red | 1 | 272.86 | 181.90 |
| 1197B20 | Disposal Can, Polyethylene, 5 Gallon, White | 1 | 282.96 | 181.90 |
| 1197B21 | Disposal Can, Polyethylene, 5 Gallon, Yellow | 1 | 282.96 | 181.90 |

3M
Virtua™

Protective Eyewear

- Three size options to fit a variety of facial profiles
- All stations are ANSI Z308, 1-2003 compliant



| Thomas No. | Description | Each | List | Sale | Case | List | Sale |
|------------|-------------------------------------|------|------|------|------|-------|-------|
| 1220W32 | Clear Temples, Clear Hard Coat Lens | 1 | 3.41 | 1.40 | 20 | 68.20 | 28.00 |
| 1186W35 | Clear Temples, Clear Hard Coat Lens | 1 | 4.02 | 2.04 | — | — | — |
| 1184L66 | Clear Temples, Clear Hard Coat Lens | 1 | 2.98 | 1.95 | — | — | — |

Honeywell

Banded Ear Plugs

- Highest attenuation among banded hearing protectors
- Soft ear pads mold quickly for a comfortable, custom fit



| Thomas No. | Description | Each | List | Sale |
|------------|-----------------------------|------|------|------|
| 8245A01 | Banded Earplugs, NRR 2 | 1 | 9.18 | 4.60 |
| 8245A03 | Replacement Ear Plugs, Pair | 1 | 5.39 | 2.61 |

MICROFLEX™
THE MOST TRUSTED NAME IN GLOVES®

SafeGrip® Latex Gloves

- Blue, textured fingers
- Ambidextrous and Powder-free



| Thomas No. | Description | Box | List | Sale | Case | List | Sale |
|------------|---------------------------|-----|-------|-------|------|--------|--------|
| 5761W21 | SafeGrip, Gloves, Small | 50 | 24.02 | 15.02 | 500 | 216.20 | 150.20 |
| 5761W26 | SafeGrip, Gloves, Medium | 50 | 24.02 | 15.02 | 500 | 216.20 | 150.20 |
| 5761W31 | SafeGrip, Gloves, Large | 50 | 24.02 | 15.02 | 500 | 216.20 | 150.20 |
| 5761W36 | SafeGrip, Gloves, X-Large | 50 | 24.02 | 15.02 | 500 | 216.20 | 150.20 |

Honeywell
NORTH

First Aid Stations

- Ideal for office, lab or manufacturing settings
- All stations are ANSI Z308, 1-2003 compliant



| Thomas No. | Description | Each | List | Sale |
|------------|---------------------------------------|------|--------|--------|
| 1233T96 | First Aid Station, 100 Person | 1 | 184.55 | 131.81 |
| 1233T97 | Lab Soaker, Light Duty Pad, 20" X 16" | 1 | 265.51 | 189.64 |

Bulk First Aid Kits

- Contains economic bulk supplies that can be readily accessed when needed
- Metal kits are made of steel, plastic kits are waterproof




| Thomas No. | Description | Each | List | Sale |
|------------|--|------|--------|-------|
| 1233T99 | Bulk First Aid Kit, 10 Person, Metal | 1 | 40.26 | 28.76 |
| 1233U00 | Bulk First Aid Kit, 25 Person, Metal | 1 | 42.30 | 30.21 |
| 1233U02 | Bulk First Aid Kit, 75 Person, Metal | 1 | 106.38 | 75.99 |
| 1233U03 | Bulk First Aid Kit, 10 Person, Plastic | 1 | 24.25 | 17.31 |
| 1233U04 | Bulk First Aid Kit, 25 Person, Plastic | 1 | 28.46 | 20.33 |
| 1233U05 | Bulk First Aid Kit, 50 Person, Plastic | 1 | 56.47 | 40.33 |
| 1233U06 | Bulk First Aid Kit, 75 Person, Plastic | 1 | 86.36 | 61.69 |

NPS CORP

Multi-Purpose Lab Soakers

- Spillfyter® Soakers provide absorption protection against spills and cushioning to minimize breakage
- Made with multiple layers of 100% recycled tissue and a plastic poly backing for diaper-like characteristics
- Heavy duty models feature durable fabric top



| Thomas No. | Description | Case | List | Sale |
|------------|---|------|--------|--------|
| 1233K17 | Light Duty Roll, 400'L x 20"W | 2 | 186.84 | 153.90 |
| 1233K18 | Light Duty Pad, 20"L x 16"W | 250 | 85.76 | 70.68 |
| 1233K19 | Heavy Duty Roll w/ Fabric Top, 200'L x 20"W | 2 | 151.36 | 124.66 |
| 1233K20 | Heavy Duty Pad w/ Fabric Top, 20"L x 16"W | 200 | 112.01 | 92.33 |

 **Kimberly-Clark**

JACKSON SAFETY* G80 NITRILE Chemical Resistant Gloves

- Provide protection against acids, caustics, bases, solvents and other aqueous solutions
- Enhanced chemical resistance, improved durability and comfort, as well as environmentally improved packaging
- Textured fingertips and palm for enhanced grip
- Now approved for food contact in the US and Canada



| Thomas No. | Description | Case | List | Sale |
|------------|--|------|--------|--------|
| 1202A89 | JACKSON SAFETY* G80 NITRILE Chemical Resistant Gloves, Size 7 - SM | 60 | 145.84 | 110.17 |
| 1202A90 | JACKSON SAFETY* G80 NITRILE Chemical Resistant Gloves, Size 8 - MD | 60 | 145.84 | 110.17 |
| 1202A91 | JACKSON SAFETY* G80 NITRILE Chemical Resistant Gloves, Size 9 - LG | 60 | 145.84 | 110.17 |
| 1202A92 | JACKSON SAFETY* G80 NITRILE Chemical Resistant Gloves, Size 10 - XL | 60 | 145.84 | 110.17 |
| 1202A93 | JACKSON SAFETY* G80 NITRILE Chemical Resistant Gloves, Size 11 - XXL | 60 | 145.84 | 110.17 |

slice™

Auto-Retractable Box Cutter

- Auto-retractable for extra safety
- High visibility green rubberized interior grip
- Includes one double-sided ceramic blade



| Thomas No. | Description | Case | List | Sale |
|------------|---------------------------------------|------|--------|--------|
| 1191X86 | Box Cutter, Auto-Retractable, Ceramic | 6 | 203.20 | 139.90 |

Manual Pen Cutter

- Manual 3-position slider button
- Includes one double-sided ceramic blade
- Stays sharp up to 10X longer than steel



| Thomas No. | Description | Case | List | Sale |
|------------|---|------|--------|--------|
| 1191X98 | Pen Cutter, 3 Position, Manual, Ceramic | 12 | 264.23 | 181.87 |

Order Now!

Inside US **800.524.0018** Outside US **856.467.2000** | thomassci.com | global@thomassci.com
Use **Promo Code 0217** to receive the featured pricing. *Valid until 12/31/17.*



Right-to-Know, Safety-Vented Wash Bottles w/GHS Labeling

- Wide mouth 53mm (2") bottle for easy refilling
- Color-coded polypropylene closures aid in quick and easy identification
- Unique valve prevents pressure build-up eliminating potentially dangerous chemical drips



| Thomas No. | Description | Pack | List | Sale | Case | List | Sale |
|------------|---|------|-------|--------------|------|--------|---------------|
| 1199R27 | Acetone, 500mL | 4 | 58.00 | 32.78 | 24 | 348.00 | 196.68 |
| 1199R28 | Dichloromethane, 500mL | 4 | 52.26 | 46.52 | — | — | — |
| 1199R29 | Deionized Water, 500mL | 4 | 58.00 | 32.78 | 24 | 348.00 | 196.68 |
| 1199R30 | Distilled Water, 500mL | 4 | 52.26 | 46.52 | — | — | — |
| 1199R31 | Ethyl Acetate, 500mL | 4 | 52.26 | 46.52 | — | — | — |
| 1199R32 | Isopropanol, 500mL | 4 | 58.00 | 32.78 | 24 | 348.00 | 196.68 |
| 1199R33 | Machine Oil, 500mL | 4 | 52.26 | 46.52 | — | — | — |
| 1199R34 | Methanol, 500mL | 4 | 58.00 | 32.78 | — | — | — |
| 1199R35 | Methyl Ethyl Ketone, Green Cap, 500mL, Non-Vented | 4 | 52.26 | 46.52 | — | — | — |
| 1199R36 | Isotonic Saline, 500mL | 4 | 52.26 | 46.52 | — | — | — |
| 1199R37 | Bleach, 500mL | 4 | 52.26 | 46.52 | — | — | — |
| 1199R38 | Toluene, Red Cap, 500mL, Non-Vented | 4 | 52.26 | 46.52 | — | — | — |
| 1199R39 | Water, 500mL | 4 | 52.26 | 46.52 | — | — | — |
| 1199R40 | Ethanol, 500mL | 4 | 58.00 | 32.78 | — | — | — |
| 1199R41 | Assortment, 500mL | 4 | 52.26 | 46.52 | — | — | — |
| 1199R42 | Acetone, 1000mL | 4 | 39.00 | 34.71 | — | — | — |
| 1199R43 | Distilled Water, 1000mL | 4 | 39.00 | 34.71 | — | — | — |
| 1199R44 | Isopropanol, 1000mL | 4 | 39.00 | 34.71 | — | — | — |
| 1199R45 | Methanol, 1000mL | 4 | 39.00 | 34.71 | — | — | — |
| 1199R46 | Bleach, 1000mL | 4 | 39.00 | 34.71 | — | — | — |
| 1199R47 | Ethanol, 1000mL | 4 | 39.00 | 34.71 | — | — | — |
| 1199R48 | Assortment, 1000mL | 4 | 77.93 | 69.37 | — | — | — |



SemperShield® Nitrile Gloves

- Extra thickness and an extended cuff
- Excellent chemical resistance
- Tested for use with chemotherapy drugs
- Manufactured without natural rubber
- Powder-Free



| Thomas No. | Description | Box | List | Sale | Case | List | Sale |
|------------|-------------------------------|-----|------|-------------|------|-------|--------------|
| 1186P92 | 5 mil, 12" Extended Cuff, S | 50 | 6.74 | 6.49 | 500 | 67.40 | 64.90 |
| 1186P93 | 5 mil, 12" Extended Cuff, M | 50 | 6.74 | 6.49 | 500 | 67.40 | 64.90 |
| 1186P94 | 5 mil, 12" Extended Cuff, L | 50 | 6.74 | 6.49 | 500 | 67.40 | 64.90 |
| 1186P95 | 5 mil, 12" Extended Cuff, XL | 50 | 6.74 | 6.49 | 500 | 67.40 | 64.90 |
| 1186P96 | 5 mil, 12" Extended Cuff, XXL | 50 | 6.74 | 6.49 | 500 | 67.40 | 64.90 |



XCEED® Powder-Free Nitrile Exam Gloves



- **Unparalleled Comfort:** In a comprehensive study XCEED® gloves markedly reduced muscle effort
- **Improved Strength:** Up to 60% stronger than leading nitrile brands and up to 90% stronger than leading latex brands
- **Leading Protection:** Setting the standard with pinhole rates that are nearly 75% below the already stringent level set by the FDA
- **Value:** Packed 250 gloves per box (230 for X-Large), XCEED® gloves save you money, time and storage space

Traditionally, glove users have had to compromise strength in order to receive a more comfortable lighter weight glove. Less comfortable gloves can increase the muscle effort required to perform repetitive tasks. Prolonged muscle effort over longer periods of time can cause hand fatigue which has been associated with more serious injuries.

| Thomas No. | Description | Box | List | Sale | Case | List | Sale |
|------------|-------------------------------------|-----|-------|--------------|------|--------|---------------|
| 1220U55 | XCEED® Powder-Free Nitrile Exam, XS | 250 | 33.29 | 19.62 | 2500 | 299.61 | 196.20 |
| 1220U56 | XCEED® Powder-Free Nitrile Exam, S | 250 | 33.29 | 19.62 | 2500 | 299.61 | 196.20 |
| 1220U57 | XCEED® Powder-Free Nitrile Exam, M | 250 | 33.29 | 19.62 | 2500 | 299.61 | 196.20 |
| 1220U58 | XCEED® Powder-Free Nitrile Exam, L | 250 | 33.29 | 19.62 | 2500 | 299.61 | 196.20 |
| 1220U59 | XCEED® Powder-Free Nitrile Exam, XL | 230 | 33.29 | 19.62 | 2300 | 299.61 | 196.20 |



Simple Green Safety Towels®

- Specially designed with the industrial cleaning strength of Crystal Simple Green infused into thermally bonded, dual textured, melt-blown, lint-free polypropylene towels
 - Safe and effective way to remove grease, oil, wet paint, carbon, graphite, tar, soot, grime, most inks and other soils from: metals finished wood, laminates, fiberglass, porcelain, tile, ceramics, concrete, terrazzo, carpets ... and other industrial surfaces
 - Pre-moistened with non-toxic and biodegradable Crystal Simple Green®
- Multi-purpose cleaning towels can handle just about any clean-up challenge.



| Thomas No. | Description | Case | List | Sale |
|------------|--------------------------------------|------|--------|---------------|
| 1200X47 | Simple Green Safety Towels, 75 Count | 6 | 167.39 | 108.94 |

Lemon Scent Simple Green® All-Purpose Cleaner

- Original non-toxic, biodegradable, concentrated formula is now available in lemon scent
- Can be diluted with water to remove grease, soap scum, coffee and juice stains, ink, lipstick and more from any washable, indoor or outdoor surface
- Lemon Scent is ideal for kitchens and baths!



| Thomas No. | Description | Case | List | Sale |
|------------|--|------|--------|---------------|
| 1200X70 | Lemon Scent Simple Green All-Purpose Cleaner, 24 oz. | 12 | 139.18 | 106.30 |
| 1200X71 | Lemon Scent Simple Green All-Purpose Cleaner, 1 gal. | 6 | 143.45 | 109.47 |



We Believe You Are Important
— How Can We Help?

Inside US 800.524.0018
Outside US 856.467.2000 • Fax 856.467.0512
PO Box 99 • Swedesboro, NJ 08085
thomassci.com • global@thomassci.com

