

Thomas BRAND

Mini DryBath

A personal sized drybath

- Compact, fits in the palm of your hand
- Exchangeable blocks, for tube sizes 0.2 to 50ml
- Digital Temperature Control
- Portable (Optional 12V car adapter)



The Thomas® Mini DryBath is a personal sized instrument that fits in the palm of your hand. Six block configurations are available for a variety of tube sizes, ranging from 0.2 to 50ml.

A large, backlit display is used to control temperature up to 80°C (minimum temperature is ambient +5°C). Each unit includes a clear plastic cover that provides enhanced temperature uniformity and ramping time. A convenient block lifter is also included, allowing block to be removed from the instrument before it has cooled.

Specifications	
Temperature Range:	Ambient +5 to 80°C
Temperature Accuracy:	+/- 0.5°C
Temperature Increments:	0.1°C
Temperature Uniformity:	+/- 0.2°C
Block Construction:	High Grade Aluminum
Dimensions:	4.4 x 5.9 x 4 in. (11.2 x 15 x 10 cm)
Weight:	1.2 lbs / 3kg
Electrical:	100 to 240V 50-60 Hz

Thomas No.	Description	Each
1222U43	Thomas® Mini DryBath, 100-240V US Plug	1
1222U44	Thomas® Mini DryBath, 100-240V EU Plug	1
1222U45	Block, 40 x 0.2ml tubes or 5 PCR strips	1
1222U46	Block, 24 x 0.5ml centrifuge tubes	1
1222U47	Block, 15 x 1.5ml centrifuge tubes	1
1222U48	Block, 15 x 1.5ml or 2.0ml tubes	1
1222U49	Block, 4 x 15ml centrifuge tubes	1
1222U50	Block, 2 x 50ml centrifuge tubes	1
1222U51	12V vehicle power adapter	1