

PlateSeal Films and Foils

Selection Guide by Application

Product Description

(Click on product name to order Free Samples)

Product Description (Click on product name to order Free Samples)	Catalog #s	Plate / Tube Type					Application						
		Raised-Rim Microplates	Flat-Top Microplates	PCR 8-Tube Strips / Strip Wells / Single Plate Row	qPCR (Real-Time)	PCR	ELISA / Incubation / General Purpose	Cell / Tissue / Bacterial Culture	Long-Term Storage	Automation Compatible	Protein Crystallography	HPLC / Chromatography	Sonication
PlateSeal 1X8 PCR SealStrips™	1150P96		*	*	*	*							
PlateSeal Clear Polyester Film	1150P94	*	*		*	*							
PlateSeal Clear Polyester Film, X-Large	1150P95		*		*	*							
PlateSeal Optically Clear Polyolefin Film, Pressure-Activated Silicone Adhesive	1150P97	*	*		*			*		*	*		
PlateSeal Optically Clear Polyolefin Film, Pressure-Activated Silicone Adhesive, Sterile	1150P98	*	*		*			*		*	*		
PlateSeal 1X8 Foil SealStrips™	1150Q04		*	*		*			*				
PlateSeal Chemical Resistant Foil	1150Q01	*	*			*		*	*				
PlateSeal Chemical Resistant Foil, Sterile	1150Q02	*	*			*		*	*				
PlateSeal QuickApply™ Dispenser with Pre-Scored Chemical Resistant Foil Roll	1150Q03	*	*			*		*	*				
PlateSeal Easy-to-Pierce Foil	1150P99		*			*		*	*				
PlateSeal Easy-to-Pierce Foil, Sterile	1150Q00		*			*		*	*				
PlateSeal 96-Well Pre-Slit Pierceable Film	1150Q09		*						*		*	*	
PlateSeal 96-Well Pre-Slit Pierceable Film, Sterile	1150Q10		*						*		*	*	
PlateSeal 96-Well Adhesive-Free Pierceable Film	1150Q11		*						*		*		
PlateSeal 1X8 Incubation SealStrips™	1150P92		*	*			*						
PlateSeal 1X8 Incubation SealStrips™, Sterile	1150P93		*	*			*						
PlateSeal Incubation Film	1150P87		*				*						
PlateSeal Incubation Film, Sterile	1150P88		*				*						
PlateSeal QuickSeal™ Incubation Film Applicator System	1150P89		*				*						
PlateSeal QuickSeal™ Incubation Film Replacement Rolls	1150P90		*				*						
PlateSeal QuickApply™ Dispenser with Pre-Scored Film Roll	1150P91		*				*						
PlateSeal Breathable Film	1150Q05		*				*						
PlateSeal Breathable Film, Sterile	1150Q06		*				*						
PlateSeal QuickSeal™ Breathable Film Applicator System, Sterile	1150Q07		*				*						
PlateSeal QuickSeal™ Breathable Film Replacement Rolls, Sterile	1150Q08		*				*						