

SPECIAL OFFER

Scale up your cell counting efficiency and throughput right from the start with the Corning® Cell Counter

Receive 4 cases of disposable slides and 2 bottles of Trypan Blue with every Corning Cell Counter purchase*.

May 15, 2022 through December 31, 2022

The Corning Cell Counter is the first cell counter that combines the best features of automated and manual cell counting. It is fast, accurate, and affordable. Now you can purchase associated consumables to scale.

Corning Cell Counter

- **Evolving** – Periodic, seamless upgrades, added features, and improved functionality.
- **Accurate** – Thanks to its cloud-based machine learning algorithm.
- **Improved resolution** – Ability to accurately count as small as 4 μm (e.g., PBMC).
- **Low-cost** – Works with common reusable glass hemocytometers. No consumables required.
- **Fast** – Thanks to its online image processing.



Corning Counting Chamber

Disposable counting chambers are composed of two 10 μL chambers. Designed to optimize throughput by enhancing speed and data accuracy in mammalian cell counting.



Corning Trypan Blue Solution

Trypan Blue is the most common stain used in cell count and viability assays. Viable cells exclude the dye, while nonviable cells absorb the dye and appear blue.



► Purchase a Corning Cell Counter (see below)...

Thomas No.	Corning Cat. No.	Description	Qty/Cs
1159B19	6749	Corning Cell Counter	1

► ...Receive at no additional cost (see below)*

Thomas No.	Corning Cat. No.	Redeem	Description
22A00D422	480202	4 cases	Disposable Counting Chamber, 0.1 mm (50/Cs)
B003L40	25-900-CI	2 cases	100 mL Trypan blue solution, 0.4% (w/v) in PBS, pH 7.5 \pm 0.5

For purchases May 15, 2022 through December 31, 2022

To redeem this promotion:

- Redeem this promotion by visiting www.corning.com/promotions
- Complete the web form, and use promotional code "PRO-CLS3015T"

Terms and Conditions: Offer valid in the U.S. and Canada only and void where prohibited by government, law, or company policy. Offer valid for end-user customers only for purchases from May 15, 2022 through December 31, 2022 or while supply lasts. Redemption information must be received by January 15, 2023. Corning has the right to substitute products if required. Please allow 4 to 8 weeks for delivery. ***Customers are limited to 10 redemptions.**

Warranty/Disclaimer: Unless otherwise specified, all products are for research use or general laboratory use only.* Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. These products are not intended to mitigate the presence of microorganisms on surfaces or in the environment, where such organisms can be deleterious to humans or the environment. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications. *For a listing of US medical devices, regulatory classifications or specific information on claims, visit www.corning.com/resources.

Corning's products are not specifically designed and tested for diagnostic testing. Many Corning products, though not specific for diagnostic testing, can be used in the workflow and preparation of the test at the customers discretion. Customers may use these products to support their claims. We cannot make any claims or statements that our products are approved for diagnostic testing either directly or indirectly. The customer is responsible for any testing, validation, and/or regulatory submissions that may be required to support the safety and efficacy of their intended application.

For a listing of trademarks, visit www.corning.com/clstrademarks. All other trademarks are the property of their respective owners.
© 2022 Corning Incorporated. All rights reserved. 5/22 PRO-CLS3015