

**Thomas**  
**BRAND**

# Mini Centrifuge

- Eight Position Microtube rotor
- Storage compartment for PCR rotor
- Near silent operation
- Starts and stops in seconds
- (8) x 0.5mL Tube adapters for Microfuge Tube rotor included with purchase
- (8) x 0.2mL Tube adapters for Microfuge Tube rotor included with purchase



The Thomas® Mini Centrifuge is a personal centrifuge that virtually fits in the palm of your hand. Ideal for quick spin downs of microtubes and PCR tubes, the Thomas® Mini Centrifuge is extremely easy to use. Simply close the lid and the rotor instantly reaches 6,000 rpm. Open the lid, and the rotor quickly decelerates to a stop for removing samples.

An increased capacity, eight-position microtube rotor is included and comes installed. A second rotor for 0.2ml PCR tubes and strips is also included and stored with the rotor removal key in a convenient compartment on the bottom of the centrifuge. When the rotor is not in use, storage in this compartment prevents loss or misplacement of the rotor.

| Thomas No.     | Description                                       | Each |
|----------------|---|------|
| <b>1222U60</b> | Thomas® Mini Centrifuge, Red lid 100-240V US Plug | 1    |
| <b>1222U61</b> | Thomas® Mini Centrifuge, Red lid 100-240V EU Plug | 1    |