

# **Technical Data**

D.C.L.S. Agar

D.C.L.S. Agar is used as a selective medium for the isolation of *Shigella* and *Salmonella* . Also useful for isolation of *Vibrio cholerae* .

## Composition\*\*

Ingredients	Gms / Litre
Proteose peptone	7.000
Beef extract	3.000
Lactose	5.000
Sucrose	5.000
Sodium citrate	10.000
Sodium thiosulphate	5.000
Sodium deoxycholate	2.500
Neutral red	0.030
Agar	12.000
Final pH ( at 25°C)	7.2±0.2

<sup>\*\*</sup>Formula adjusted, standardized to suit performance parameters

#### **Directions**

Suspend 49.53 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. DO NOT AUTOCLAVE. Cool to 50°C and pour about 20 ml of medium into sterile Petri plates and allow to dry for about two hours with covers partially removed.

# **Principle And Interpretation**

Salmonella infection leads to salmonellosis, which ranges clinically from self-limited gastroenteritis (diarrhea, abdominal cramps and fever) to enteric fevers (including typhoid fever). Shigella species cause classical bacillary dysentery characterized by severe cramping abdominal pain and diarrhea with blood and mucus.

D.C.L.S. Agar is a modification of Deoxycholate Citrate Agar of Leifson (1). It is a slightly selective and differential medium, which incorporates sucrose as an additional fermentable carbohydrate to differentiate lactose negative sucrose positive coliforms from *Salmonella* species The addition of sucrose to this medium increases its usefulness because non-pathogenic sucrose fermenting organisms like *Proteus*, *Enterobacter*, *Klebsiella* form red colonies. D.C.L.S. Agar is a moderately selective culture medium which also supports the growth of *Vibrio* species.

## **Quality Control**

### **Appearance**

Light yellow to pink homogeneous free flowing powder

#### Gelling

Firm, comparable with 1.2% Agar gel.

#### Colour and Clarity of prepared medium

Reddish orange coloured, clear to slightly opalescent gel forms in Petri plates

#### Reaction

Reaction of 4.95% w/v aqueous solution at 25°C. pH: 7.2±0.2

#### рH

7.00-7.40

# **Cultural Response**

M160: Cultural characteristics observed after an incubation at 35-37°C for 18-48 hours.

Organism	Inoculum	Growth	Recovery	Colour of
	(CFU)			colony

HiMedia Laboratories Technical Data

Enterococcus faecalis ATC	$C >= 10^3$	inhibited	0%	
Escherichia coli ATCC 25922	50-100	none-poor	<=10%	red
Proteus vulgaris ATCC 13315	50-100	luxuriant	>=50%	red
Salmonella Typhimurium ATCC 14028	50-100	luxuriant	>=50%	colourless- slightly pink
Shigella flexneri ATCC 12022	50-100	fair-good	30-40%	colourless- slightly pink

## **Storage and Shelf Life**

Store below 30°C in tightly closed container and the prepared medium at 2 - 8°C. Use before expiry date on the label.

#### Reference

1. Leifson E., 1935, J. Pathol. Bacteriol., 40:581.

Revision: 1 / 2011

# (

#### Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.