

# Weighing Supplies

*A variety of styles to suit your weighing requirements*

## Aluminum Moisture Pans

*Ideal for loss and drying determinations*

- Withstand temperatures up to 677°C (1,250°F)
- Can be used as evaporating dishes or dust covers
- Offered with a flat bottom, slightly crimped side walls and a small tab handle, or with smooth walls and no tab handle



9885D01



9885D19

Thomas No.	Size (mm)	Diam. (In.)	Approx. Cap.	Depth (mm/in.)	Weight	Pack	Case
<b>Aluminum Weigh Dishes with Crimped Side Walls and Tab Handle</b>							
<b>9885D01</b>	43	1.75	20 mL	7 mm / 5/16 in.	1.0 g	100	1000
<b>9885D07</b>	57	2.25	50 mL	16 mm / 5/8 in.	1.5 g	100	1000
<b>9885D13</b>	70	2.75	70 mL	17 mm / 11/16 in.	2.0 g	100	1000
<b>Aluminum Weigh Dishes with Smooth Walls and No Tab Handle</b>							
<b>9885D19</b>	57	2.25	50 mL	16 mm / 5/8 in.	1.5 g	100	1000
<b>9885D28</b>	70	2.75	70 mL	17 mm / 11/16 in.	2.0 g	100	1000

## Pre-Cut Aluminum Squares *Eliminate time spent cutting aluminum foil*

- Multiple uses: Cap flasks, beakers, bottles and other labware
- Act as a barrier or insulator
- Easy to separate: Squares are interleaved with paper
- High resistance to fats, petroleum, greases and organic solutions non-toxic, non-magnetic, non-sparking and corrosion resistant; good electrical and thermal conductivity

Thomas No.	Size (In.)	Thickness (In.)	No. of Squares*
<b>1086J24</b>	4 x 4	0.001	2200
<b>1086J26</b>	6 x 6	0.001	1100
<b>1086J28</b>	8 x 8	0.0015	430
<b>1086J30</b>	6 x 6	0.0015	192
<b>1203N63</b>	12 x 12	0.0015	1000

\*Furnished in 5 lb. packages; Count is approximate



*Continued on reverse side*

## Weighing Dishes

- Available in Hexagonal and Square
- Flat bottom
- Perfect for weighing solids and liquids
- Made from flexible polystyrene
- Provides excellent protection against spills
- Prevents hazards and sample loss



Thomas No.	Size	Dimensions (In.)	Pack	Case
<b>Square</b>				
<b>9885D37</b>	39 mm	1-1/2 x 1-1/2 x 3/8	500	6000
<b>9885D40</b>	77 mm	3 x 3 x 1	500	2000
<b>9885D43</b>	127 mm	5 x 5 x 1	500	1000



Thomas No.	Dimensions (In.)	Box	Case
<b>Hexagonal</b>			
<b>9885D55</b>	1-1/2 x 1	500	4000
<b>9885D58</b>	2-1/2 x 2	500	3000
<b>9885D60</b>	4 x 3	500	2000
<b>9885D63</b>	5 x 3-1/2	500	2000




## Weighing Papers

- Economical
- Smooth, strong and non-absorbent with glossy surface
- Especially suited for weighing powders

Thomas No.	Dimensions (In.)	Box
<b>9885G40</b>	3 x 3	500
<b>9885G50</b>	4 x 4	500
<b>9885G60</b>	6 x 6	500

## Weighing Canoes



- Aerospace-grade, anti-static plastic
- Natural opaque white color
- Economical, stable, easy-to-handle design avoids "squeezing" of square weigh boats
- Plastic provides reduced adhesion without oily coatings

Thomas No.	Dimensions (In.)	Pack	Case
<b>9885E07</b>	2-1/4 x 1 x 1/4	250	3000
<b>9885E09</b>	5 x 3-1/4 x 1	250	1500
<b>9885E11</b>	7 x 4-1/4 x 1-1/8	250	1000