

Clip COVID Rapid Antigen Test & Analyzer - Nasal Swab - Point-of-Care - 200 Tests

Quantity: 200 Test Kits (25 ea x 2) and Analyzer(s) Placement

Description: FDA-authorized easy to use 30-minute rapid antigen test for SARS-CoV-2. Get a readout on an iPhone-based analyzer with 100% negative agreement (specificity) and 96.9% positive agreement (sensitivity) compared to RT-PCR. An objective positive/negative readout. Only 90 seconds of hands on time per test. To learn more, watch the [3 Minute Demonstration](#). Discounted per-test pricing, down to \$23 per test, available for higher quantities (click Show More below).

Collection Method/Analysis Location: CLIA Waived/POC

Sample: Anterior nasal swab

Diagnostic Target: Antigen

Time to Results: 30 minutes

What's in each box:

Clip COVID Test Kit (MDL002): Cartridges (25), Extraction Tubes (25), pre-filled with assay reagents, Dropper Tips (25), Sterile Nasal Swabs (25), Positive Control Swab (1), Negative Control Swab (1), Package Insert (1), Quick Reference Instructions (1)

Clip Analyzer (MDL001): Clip Analyzer (1), Charging Cord (1), Power Adapter (1), User Manual (1), Tube Holder (1)

Additional Resources:

[Instructions for Use/Package Insert](#)

[Quick Reference Instructions](#)

[Clip Analyzer User Manual](#)

[FDA Letter of Authorization](#)

[Fact Sheet for Healthcare Providers](#)

[Fact Sheet for Patients](#)

Reimbursement under CPT code 87426 (Medicare rate \$45.23)

Discounted Pricing (includes Clip Analyzer Placement for New Customers)

- \$31 per test for 200 Tests (Quantity = 1 below)
- \$29 per test for 400 to 600 Tests (Quantity = 2 to 3 below)
- \$28 per test for 800 to 1,000 Tests (Quantity = 4 to 5 below)
- \$27 per test for 1,200 to 1,400 Tests (Quantity = 6 to 7 below)
- \$26 per test for 1,600 to 1,800 Tests (Quantity = 8 to 9 below)
- \$25 per test for 2,000 to 2,800 Tests (Quantity = 10 to 14 below)
- \$24 per test for 3,000 to 4,800 Tests (Quantity = 15 to 24 below)
- \$23 per test for 5,000 to 10,000 Tests (Quantity = 25 to 50 below)

New Customers will need to enter into an Analyzer Placement Agreement with Luminostics before their first order can be fulfilled. 1 Clip Analyzer will be placed at no charge with your first order, and additional Clip Analyzers are available for a nominal additional cost.

Per the CDC: For the best results, collect and handle the sample according to the manufacturer's instructions. If you don't collect the sample as directed, your test results may be incorrect. Once collected, use the sample as described in the manufacturer's instructions.

This test has not been FDA cleared or approved, but has been authorized by FDA under an EUA for use by authorized laboratories; This test has been authorized only for the detection of proteins from SARS-CoV-2, not for any other viruses or pathogens; and This test is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostics for detection and/or diagnosis of COVID-19 under Section 564(b)(1) of the Federal Food, Drug and Cosmetic Act, 21 U.S.C. 360bbb-3(b)(1), unless the declaration is terminated or authorization is revoked sooner.